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GB Nelaton Catheter - Instruction for Use

DESCRIPTION
The Nelaton Catheter is designed for transient use and temporary drainage of urine. Its flexible tube is placed in the body to drain or collect urine from the bladder. It is made from medical grade PVC, consists of shaft and connector. The Nelaton Catheter is used in various departments, such as department of urology, internal medi- cine, surgery, obstetrics, and gynecology. It is used for patients suffering from moving difficulty or being completely bed-ridden.

FEATURES
- Color coding for size identification
- Distal rounded end for non-traumatic introduction
- Available with heat punching lateral eyes for efficient drainage and smooth painless insertion
- Available with male, female and tiemann tip.

INTENDED USE / INDICATIONS
Nelaton catheter is inserted into the bladder through urethra and are indicated for intermittent urine drainage.

CONTRAINDICATIONS
- Acute urethritis
- Acute prostatitis
- Acute epididymitis

SUITABLE PATIENTS:
- Infant: 6Fr ;
- Pediatric: 8Fr, 10Fr and 12Fr ;
- Adult: 12Fr, 14Fr, 16Fr, 18Fr, 20Fr, 22Fr and 24Fr

ADVERSE REACTIONS
- Catheter draining little or no urine despite enough fluid intake
- Leakage of large amounts of urine around the catheter
- Urinary tract or kidney infections
- Urethritis
- Urethral Bleeding or Injury

DIRECTIONS FOR USE
- Always wash hands before use.
- Using aseptic technique, remove the catheter from its pouch and place it on a sterile field.
- Peel open the package and take out the catheter.
- Lubricate the catheter with water-base lubricant.
- Put on sterile gloves, hold the connector of the catheter with the dominant hand and insert the catheter slowly and gently into the urethra until urine begins to flow. Direct the connector to a suitable container or toilet.
- When the urine stops flowing, gently press the bladder position to make sure all urine is drained out. Gently withdraw the catheter. Stop if the urine starts to flow again.
- Clean the urethral orifice and the surrounding skin.
- Put the catheter back into the package and dispose in accordance with accepted medical practice and applicable regulations. Do not flush it down to the toilet.
- Wash hands.

WARNINGS / PRECAUTIONS
- Read all warnings and instructions before use. Improper use can result in serious or fatal illness or injury.
- Strictly follow the clinical professional's instruction on interval time of urine drainage.
- This device can be operated by medical care person or patients who are trained on the drainage procedure. If operated by patient, any question during operation or use, consult medical care professionals.
- Do not use if the sterile packaging is damaged or unintentionally opened before use.
- Do not use if it's damaged or irregularly shaped.
- Do not use after the expiry date.

- For single use only.
- Choose proper size. Small catheter is recommended for those suspected of urinary tract stenosis.
- Insert the urethrain a gentle manner, notto damage the urethra mucosa. Do not insert too deep, especially avoid twitching the catheter.
- Drain the urine slowly for the patient with overfilling bladder or in extremely weak to avoid bleeding or syncope caused by excessive decompression of the bladder.
- Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in injury, illness or death of the patient. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Dispose of product and packaging in accordance with hospital administrative and/or local govern- ment policy.
- Any serious incident that has occurred in relation to the device should be reported to the manu- facturer and the local competent authority of the user's place.

STORAGE CONDITIONS
- Store in a dark, cool and dry place
- Protect product from moisture and excessive heat.
- Avoid prolonged exposure to ultraviolet, sunlight and fluorescent light.
- Store in manner preventing crushing.

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Well Lead Medical Co., Ltd. C-4 Jinhu Industrial Estate Hualong 511434 Panyu Guangzhou People's Republic of China
Shanghai International Holding Corp. GmbH Eiffestraße 80 20537 Hamburg Germany
QS Engineering AG Erlenstrasse 31 4106 Therwil Switzerland