FILTER TECHNOLOGY

Instructions For Use; Mode d'emploi; Gebrauchsanweisung; Istruzioni per l'uso; Instruções de uso; Gebruiksaanwijzing; Instruções de Utilização; Bruksanvisning; Bruksanvisning; Bruksanvisning; Käyttöohjeet; Instrukcja użytkowania; Иструкции за употреба; Návod k použití; Kasutusjuhend; Οδηγίες χρήσης; Használati útmutató; Treoracha Úsáide; Lietošanas instrukcija; Naudojimo instrukcijos; Istruzioni per l'uso; Instruções de utilização; Návod na použitie; Navodila za uporabo.

CE 2797

Table with 6 columns: Classification, UMDNS code, Packaging code, UK, FRANCE, GERMANY, SPAIN, NETHERLANDS, USA. It details regulatory compliance for various regions.

Table with 16 columns: REF, 4020/01XXXX ECO Machine Filter, 4020/03XXXX Single Walled Machine Filter, 4020/06XXXX ECO Machine Filter, 4220/01XXXX Eco maxipleat, 4220/04XXXX ECO Filter, 4244/01XXXX ECO maxi-pleat, 4244/01DBKXXXX Pleated filter and collapsible tube, 4244/01DEKXXXX Maxi pleat with connectors, 4244/02XXXX ECO Maxi Pleat, 4244/04XXXX ECO maxi-pleat, 4244/06XXXX Eco maxi-pleat, 4244/700XXXX MAXI Maxi Angled Pleat, 4244/701XXXX ECO Maxi Pleat, 4244/702XXXX MAXI Maxi Pleat, 8444/01XXXX Maxi pleat/1, 8444/27XXXX Maxi pleat. Each cell contains detailed technical specifications in multiple languages.

Table with 16 columns: STERILE EO. Each cell contains sterilization details in multiple languages, including methods like ethylene oxide and sterilization conditions.

Table with 16 columns: NON STERILE. Each cell contains non-sterile details in multiple languages, including material composition and handling instructions.

Ambient Conditions: Conditions ambiantes : Umgebungsbedingungen : Condizioni ambientali : Condiciones ambientales : Omgivingscondities : Condições ambientais : Omgivningsförhållanden : Omgivelsesbetingelser : Krav til omgivelse : Ympäristöolosuhteet : Warunki otoczenia : Услови на околната среда : Увjeti okoline : Okolní podmínky : Keskonnattingimused : Συνθήκες περιβάλλοντος : Környezeti feltételek : Coinníollacha Comhthimpeallacha : Arkáteljele apstájkai : Aplinkos sąlygos : Kundizjonijet Ambientali : Condizioni ambientali : Okolité podmínky : Okolje.

Table with 16 columns: Temperature during storage. Each cell contains storage temperature ranges in multiple languages, such as 5 to 40 °C (41-104 °F).

Performance Data: Características ventilatorias : Leistungsdaten : Dati di rendimento : Características de rendimiento : Karakteristieke waarden : Dados de desempenho : Prestandauppgifter : Ydelse : Kapasitetsdata : Suorituskykytiedot : Parametry sprężyt : Дани за ефективността : Podaci o izvedbi : Údaje o výkonu : Toimivusandmed : Δεδομένα απόδοσης : Teljesítmény adat : Sonrai Feidhmíochta : Veiktspējas dati : Pasirodymo data : Dejta tal-Prestazzjoni : Date de performanță : Údaje o výkonnosti : Podatki o uspešnosti.

Table with 16 columns: Internal Volume. Each cell contains internal volume values in multiple languages, ranging from 60ml to 74ml.

Table with 16 columns: Resistance. Each cell contains resistance values in multiple languages, such as Max. 194.7Pa @ 30L/min.

Table with 16 columns: Filter Efficiency. Each cell contains filter efficiency percentages in multiple languages, such as Min. 99.9%.

Table with 16 columns: Bacterial Efficiency. Each cell contains bacterial efficiency percentages in multiple languages, such as Min. 99.999%.

Table with 16 columns: Viral Efficiency. Each cell contains viral efficiency percentages in multiple languages, such as Min. 99.999%.

Table with 16 columns: Gas leakage. Each cell contains gas leakage values in multiple languages, such as Max. 0.0ml/min @ 7kPa.

Table with 16 columns: Tidal Volume range. Each cell contains tidal volume ranges in multiple languages, such as 150 - 1500ml.

Material Data: Matériau: Material. Dati materiale: Material; Material; Material; Materiale; Materiale; Materiali; Material; Дани за материала: Podaci o materijalu: Materialové údaje; Materialised andmed; Δεδομένα υλικού; Anyagi adatok; Sonrai Ábhar; Materialle dati; Esminial duomenys; Dejta Materjali; Date materiale; Materialové údaje; Materialni podatki.

Table with 16 columns: Housing. Each cell contains housing material information in multiple languages, such as PP, PP + SBC, Glass Microfibre Media.

English Instructions for Use

PLAID FILTERS.

Intended use: Filters protect the patients airways effectively from exogenous microbial loads, thus reducing the risk of extrinsic colonisation and infection. Used to help reduce cross contamination between patient and machine. The Luer port connector is used for monitoring respiratory and/or anaesthesia gases.

Bidirectional Filters.

Filters can be supplied with accessories, such as catheter mounts.

Intended use of catheter mounts: Catheter mounts are used in conjunction with filters in ventilator patient circuits. For performance data, refer to the Instruction For Use for associated filter.

Cautions

CAUTION: A CAUTION statement provides important information about a potentially hazardous situation which, if not avoided, may result in harm to the patient.

Any use of the medical device requires full understanding and strict observation of all portions of these Instructions for Use. The medical device may only be used for the purpose specified under "Intended Use". Observe all CAUTION statements throughout this manual and all statements on medical device labels. The manufacturer disclaims any liability for patient injury caused by noncompliance with these statements.

Following installation of the filter, check that all connections are secure.

The medical device IFU is not available for each individual filter. Only one copy of the Instructions for Use is included in the distribution package and should therefore be kept in an accessible location for users.

MD

Medical device is only to be used by healthcare professionals.


User has the responsibility to report any serious incident that has occurred in relation to the medical device, to the manufacturer.

EC REP and the competent authority of the Member State in which the user and/or the patient is established.

To avoid contamination and soiling, the product should remain packaged until ready to be used.

Do not use the product if the packaging and/or product is damaged.

If product was kept outside the stated storage temperature, product functionality cannot be guaranteed. Storage Temperature: +5°C to 40°C.

Store the product at room temperature, in dry place, free of dust and not exposed to direct sunlight.

Handle filter with care. Abusive handling can damage the filter.

For Your Safety and that of Your Patients, follow the Instruction for Use of the machine!

Before installation, check that all system components are free of obstructions and of foreign bodies. Otherwise, use is limited, or faulty operation is possible.

Inappropriate connection to other devices can result in a patient hazard.

Medical device products have been designed, tested and manufactured exclusively for disposable use and for a period of use not exceeding 24 hours following the first opening of the primary packaging.

Do not use after expiry date stated on primary packaging.

Manufacturing date Expiry date

Monitor filter performance during the use. If any changes to filter performance or appearance are observed, change the filter to new one.

Always keep spare filter ready to exchange with the one in use, if any of the below are observed:

- Increase in resistance;
- Filter is visibly soiled or damaged or if it rattles;
- Excessive condensation;
- Associated alarm.

If the product features a Luer Port, it is used exclusively for gas monitoring. The introduction of gases or liquids via the Luer Port connection can cause serious injury to the patient.

Do not sterilize.

The products must not be reused, cleaned or sterilised / not to be sterilised after use in order to re-use. Exceptions from this rule will be stated within the Instruction For Use for the particular products code if applicable. Reuse, cleaning or sterilization may lead to failure of the products, risk infection and injury to the patient.

This product is for single patient use only.

This product is latex and PHT free.

Batch code symbol.

Prescription only.

Following use, the products must be disposed off in accordance with the local hospital hygiene and waste disposal regulations for hazardous infectious clinical waste.

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