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
2022-03

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LIGASANO® grün unsteril

Instructions for use • Mode d'emploi • Istruzioni per l'uso • Instrucciones para la aplicación • Bruksanvisning

## LIGASANO®

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# LIGASANO<sup>®</sup> green unsterile Instructions for Use

**Please read these instructions carefully before use!**

## General Information:

LIGASANO<sup>®</sup> green is an elastic, expanded polyurethane (PUR), which is almost completely pervious to air and liquids - imagine a three-dimensional sieve. It has a coarse texture and rough surface. LIGASANO<sup>®</sup> green unsterile is used for concomitant wound treatment and care, if a strength and elastic material with high permeability to air and liquid is required.

**LIGASANO<sup>®</sup> green may display a yellowish discolouration, however this does not affect the quality.**

## Intended Use:

In **non-sterile application** as a **secondary dressing (without direct wound contact)** in wound treatment or in **concomitant wound treatment**, e.g. for pressure relief, or **preventively** as **pressure and friction protection**.

## Indications:

Concomitant wound therapy in case of acute and chronic wounds (e.g. decubital ulcers, ulcus cruris, diabetic foot ulcers). Also post-surgical wounds and post-surgical wound healing impairments and thermal wounds, such as burns, scalds, frostbite.

Preventive use when the skin is still intact, e.g. pressure ulcer prophylaxis through active positioning on LIGASANO<sup>®</sup> green unsterile. Prophylaxis of maceration.

## Adverse Effects, Side Effects, Residual Risks:

There are no particular risks when using LIGASANO<sup>®</sup> green unsterile for pressure relief and hygiene, if the general nursing rules and our instructions are noted.

LIGASANO<sup>®</sup> green unsterile is a very simple product and its effects are due to its structure. It does not contain nor emit any active compounds. With the correct application of LIGASANO<sup>®</sup> green there are no known negative side effects, nor incompatibility or interactions with medicine.

**Please report any incidents that have occurred in connection with our dressing material LIGASANO<sup>®</sup> green unsterile to us or to the competent authority of the member state in which the incident occurred.**

## Contraindications:

Direct contact to skin, direct contact to wounds.

## Patient Target Group:

Children (> 2 - 16 years of age) and adults (> 16 years of age) with or without wounds.

## User/Practitioner:

**Application must always be done by healthcare professionals.** The decision about the particular kind of application of LIGASANO® green unsterile is determined by healthcare professionals and is adapted and appropriated in every single case.

## Capability Characteristics and Clinical Benefit:

With the unsterile wound dressing LIGASANO® green you can prevent the risk of the development of pressure ulcer and moisture-caused wounds.

The expected clinical benefit present itself in a very effective pressure relief. Both in already existing wounds and preventive. An additional benefit is the moisture-regulating effect.

## Application of LIGASANO® green unsterile Dressing for Pressure Relief and Prevention:

LIGASANO® green unsterile is suitable for single use only. The product must not be used after the expiration date and is not suitable for reconditioning.

LIGASANO® green unsterile must not continuously remain on the body. Recommendation: dressing change after 1-4 days, depending upon the indication. Any applications in contact, connection or combination with any additional preparations, drugs, solutions, ointments etc., are not guaranteed by us. Applications in combination with additional mechanical, electrical or electronical apparatuses or aids are not guaranteed by us.

LIGASANO® green unsterile has cushion qualities, as known from elastic foams, but without their property of thermal insulation, which is unfavourable in nursing care. The special properties of LIGASANO® green unsterile mean that parts of the body can breathe, where dampness or suppuration would otherwise accumulate as a result of contact with other objects (e.g. flat bed surfaces, seats or back-rests etc.). As a result of gravity, large amounts of fluid can run downwards, and dampness (e.g. sweat) can evaporate. Skin maceration can be effectively prevented with this product.

At this application of LIGASANO® green unsterile the contact with the skin is avoided, because the rough surface of the material is mostly felt as unpleasant and perturbing, and for preventing lesions of the skin, which could effect by friction. In addition, without any heat insulating layer, in normal temperated ambiance, it can come to hypothermia at the affected parts of the body. For an intermediate layer between skin and LIGASANO® green unsterile are exemplary suitable: stockings, clothing, bed sheet and especially a layer of LIGASANO® white with a thickness of 0.5 or 1 cm. It is important, that the intermediate layer is likewise permeable to air.

To obtain a good pressure relief by pressure distribution, you can cut LIGASANO® green unsterile to fit with a pair of scissors.

## Sterilisation on-site and subsequent application of the thus sterilised products for wound treatment:

The non-sterile product has to be sterilised at 134 °C for five minutes with moist heat according to a validated procedure. These standardised parameters meet the requirements for adequate product safety and economy, and there are no measurable effects on the properties of LIGASANO® green.

DE

GB

FR

IT

ES

RO

SE

NL

RUS

PL

CZ

LT

GR

HU

TR

When sterilising unsterile LIGASANO® products on site (e.g. in the clinic, practice, etc.), the legal regulations and risks must be observed. In particular, we point out that if the contamination is too high (> 200 CFU), the sterilised product could be contaminated with an excessive, risk-increasing amount of pathogenic germs.

**In case of sterilisation or resterilisation you act on your own responsibility. Please also note that the indication, side effects, contraindication, performance characteristics and use change when you apply the on-site sterilised product for wound treatment.**

### **Intended Use of the Products Sterilised On-site:**

Promotion of wound healing in sterile application by cleaning, covering or filling wounds.

### **Indications of the Products Sterilised On-site:**

Acute and chronic wounds (e.g. decubital ulcers, ulcer cruris, diabetic foot ulcers). Post-surgical wounds and post-surgical wound healing impairments. Thermal wounds, such as burns, scalds, frostbite.

### **Adverse Effects, Side Effects, Residual Risks of the Products Sterilised On-site:**

**Wound cleaning:** mechanical debridement may cause pain, even with LIGASANO® green sterile.

**Wound treatment:** because of its cell structure a vascularisation (e.g. of granulation tissue) into the material is possible. You can prevent this with a timely dressing change, after 1-4 days using LIGASANO® green sterile. On passive wounds or parts of the body that receive poor blood supply, pain may occur after application. This occurs when the blood circulation is stimulated by the mechanical stimulus of the material to such an extent, that previously diminished or reduced pain returns. Initially the pain can be intense, however it is usually regulated after a few hours to days. If this effect is not desired, you should forgo its application. In case of an "open wound care" without additional cover you have to take into consideration that there exists barely protection against the ambient atmosphere. Too much pressure due to the material's prior tension, external influence, or suction may lead to compression of vessels and therefore to pressure ulcers.

If hypergranulation occurs, please treat it appropriately and switch to other dressing material groups, if necessary.

LIGASANO® green is a very simple product and its effects are due to its structure. It does not contain nor emit any active compounds. With the correct application of LIGASANO® green there are no known negative side effects, nor incompatibility or interactions with medicine. Typical side effects may include reddening of the skin upon contact and initial "itching". More intense side effects may include return of sensitivity due to stimulation of local blood circulation. Its effects are only physical and end almost immediately after discontinuation of application. As with all wound dressings, minor skin reactions such as maceration, erythema, secondary infected dermatitis, erysipelas as well as hypersensitivity reactions and pain at dressing change may occur.

**Please report any incidents that have occurred in connection with our dressing material LIGASANO® green sterile to us or to the competent authority of the member state in which the incident occurred.**

### **Contraindications of the Products Sterilised On-site:**

Tumour wounds and untreated osteomyelitis. Not suitable for contact with organs, contact to exposed blood vessels or anastomotic areas with the danger of damage, or contact with nerves.

## Patient Target Group of the Products Sterilised On-site:

Children (> 2 -16 years of age) and adults (> 16 years of age) with acute, chronic, post-surgical wounds/wound healing disturbances and thermal wounds.

## User/Practitioner of the Products Sterilised On-site:

**Application must always be done by healthcare professionals.** The decision about the particular kind of application of LIGASANO® green sterile is determined by healthcare professionals and is adapted and appropriated in every single case.

## Capability Characteristics and Clinical Benefit of the Products Sterilised On-site:

With the sterile wound dressing LIGASANO® green you can **clean, cover or pack wounds**. The product can be used both for **mechanical debridement in all phases of wound healing** and for **wound treatment during the phase of exudation and granulation**. The expected clinical benefit present itself in a very effective wound cleaning and wound bed conditioning.

## Application of On-site Sterilised LIGASANO® green for Wound Treatment:

LIGASANO® green sterile is suitable for single use only. Sterility is ensured only by intact package. The product must not be used after the expiration date and is not suitable for reconditioning or resterilisation due to its potential for contamination.

LIGASANO® green sterile must not continuously remain in the wound or on the body. The application may be repeated over a period up to 30 days, always with a new device. Dressing change after 1-4 days, depending upon the indication. Please consider that LIGASANO® green sterile must not have direct contact to organs, exposed blood vessels, anastomotic areas or nerves. It may have direct contact to mucosa.

Any applications in contact, connection or combination with any additional preparations, drugs, solutions, ointments etc., are not guaranteed by us. Applications in combination with additional mechanical, electrical or electrical apparatus or aids are not guaranteed by us.

Because of its open pores, LIGASANO® green shows no suction effect. If you use it as a wound dressing, you need to combine it with a suction source or an absorbent wound pad (e. g. LIGASANO® white sterile, superabsorbents, foil and Redon bottle, NPWT system).

### ► Wound Cleaning / Mechanical Debridement:

LIGASANO® green sterile is suited for mechanical cleaning of nearly every kind of wound (see contraindications) and during all wound phases. It is highly effective in removing very persistent, viscous, and/or moist necrotic coatings. It may collect plenty of detritus. The wound and the edge of the wound can be treated with LIGASANO® green sterile. For this treatment LIGASANO® green sterile remains in the wound only for a short time to clean it. You can treat the wound with either the twisting or the wiping technique. Please ensure that the wiping direction must always be from inside to outside.

Important: please note that LIGASANO® used for wound cleaning must not remain in the wound. If the proceeding wound dressing is to be made with LIGASANO®, a new wound pad must be applied.

### ► Wound Treatment:

LIGASANO® green sterile is usually applied directly onto or into the wound without any further steps or products (e.g. ointments or irrigation).

Additional measures may negatively affect the desired result. By the described open structure of LIGASANO® green wound irrigations are possible in spite of the filled wound. Because of its rough surface the decision about the particular kind of application of LIGASANO® green is determined by healthcare professionals and is adapted and appropriated in every single case.

**Dressing composition:** if the wound dressing composition consists of several sheets or pieces, please count the number of used dressing pieces and note it in the wound documentation, to ensure the removal of all applied pieces. During application and dressing change friction of the wound, wound bed or skin should be avoided as far as possible to prevent causing any or additional lesions.

**Dressing change frequency:** one time every 1-4 days.

**Dressing technique:** each case should be evaluated and adapted accordingly. Ensure that during dressing change, the wound is left uncovered for the shortest possible amount of time in order to avoid cooling down and to keep risk of re-infection low.

### **Storage, Transport, Shelf Life:**

Store LIGASANO® green unsterile in a dry area that is protected against UV Radiation. The original packaging offers protection against humidity and UV radiation due to the filtering plastic film, and therefore storage is only permitted in the original packaging. The product may display a yellowish discolouration, especially if improperly stored, however this does not affect the quality. Special transport conditions (e.g. cooling) are not necessary. The product is stable for 60 months after the date of manufacture. The date of manufacture and expiration can be found on the product label of the original packaging.

### **Disposal / Environmental Relevance:**

In addition to its accredited quality (QM system for medical devices), LIGAMED® medical Produkte GmbH also documents environmental awareness by validation according to EC regulation and its participation in „Umweltpakt Bayern“ (Bavarian Environmental Understanding). We will send you our current certificates upon request. All LIGASANO® products are single-substance-materials made of polyurethane, consisting of approx. 97% air, and can be disposed of with ordinary household waste. Disposal of the packaging is internationally regulated in different ways. In Germany it is regulated by a licensed company.

## Symbols Used for Labelling

	Item number		Manufacturer
	Lot number		Production date
	Use by: year / month		Distributor
	Medical device		Unsterile
	Unique device identifier		Disposable product
	Consult instructions for use		Protect from sun and UV radiation
	Do not use if the packaging is damaged		Protect from moisture

### Date of Revision:

2022 / 03

### Information about the Manufacturer:

LIGAMED® medical Produkte GmbH • Pfannenstielstr. 2 • 90556 Cadolzburg/Germany • Tel. +49 9103 2046 • E-Mail info@ligamed.de

An electronic version of this instruction for use is available on our website: [www.ligasano.com/downloads](http://www.ligasano.com/downloads)

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A.Pf.E.L.®,  are Trade Marks of LIGAMED® medical Produkte GmbH, Germany
Copyright by LIGAMED® medical Produkte GmbH
**LIGAMED® medical Produkte GmbH** • Pfannenstielstr. 2 • D - 90556 Cadolzburg