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LIGASANO® white unsterile Instructions for Use

Please read these instructions carefully before use!

General Information:

LIGASANO® white is a therapeutically effective PUR foam material with a wide range of applications. The therapeutic uses are based on three underlying physical principles:

Mechanical stimulus: Stimulates local blood circulation on the skin, thereby improving the supply of nutrients and oxygen in the treated skin area and acts as a preventative on intact skin. Stimulates normal bodily functions on a local level, especially in the case of inactive patients or, where these functions are reduced. The mechanical stimulus lasts up to three days, by which time LIGASANO® white should then have been changed.

Low-pressure tension: Low-pressure adaption to contours and minimised pressure on skin and body.

Specific suction effect: Absorbs excess fluid without a drying-out effect. Averts skin maceration and is skin-friendly.

LIGASANO® white may display a yellowish discolouration, however this does not affect the quality.

Intended Use:

In **non-sterile application** as a **secondary dressing (without direct wound contact)** in wound treatment or in **concomitant wound treatment**, e.g. for pressure relief, or **preventively** as **pressure and friction protection**.

Indications:

Concomitant wound therapy in case of acute and chronic wounds (e.g. decubital ulcers, ulcerus cruris, diabetic foot ulcers). Also post-surgical wounds and post-surgical wound healing impairments and thermal wounds, such as burns, scalds, frostbite.

Preventive use when the skin is still intact, e.g. pressure ulcer prophylaxis through active positioning on LIGASANO® white unsterile. Prophylaxis of maceration and intertrigo.

Adverse Effects, Side Effects, Residual Risks:

There are no particular risks when using LIGASANO® white unsterile for pressure relief and hygiene, if the general nursing rules and our instructions are noted.

LIGASANO® white unsterile is a very simple product and its effects are due to its structure. It does not contain nor emit any active compounds. With the correct application of LIGASANO® white there are no known negative side effects, nor incompatibility or interactions with medicine. Typical side effects may include reddening of the skin upon contact and initial "itching". Its effects are only physical and end almost immediately after discontinuation of application.

Please report any incidents that have occurred in connection with our dressing material LIGASANO® white unsterile to us or to the competent authority of the member state in which the incident occurred.

Please report any incidents that have occurred in connection with our dressing material LIGASANO® white unsterile to us or to the competent authority of the member state in which the incident occurred.

Contraindications:

Very sensible skin (e.g. parchment skin), direct contact to wounds.

Patient Target Group:

Children (> 2 - 16 years of age) and adults (> 16 years of age) with or without wounds.

User/Practitioner:

Application must always be done by healthcare professionals. The decision about the particular kind of application of LIGASANO® white unsterile is determined by healthcare professionals and is adapted and appropriated in every single case.

Capability Characteristics and Clinical Benefit:

With the unsterile wound dressing LIGASANO® white you can prevent the risk of the development of pressure ulcer and moisture-caused wounds. It is able to prevent the development of intertrigo and partly fungal infections, too.

The expected clinical benefit present itself in a very effective pressure relief. Both in already existing wounds and preventive. An additional benefit is the moisture-regulating effect.

Application of LIGASANO® white unsterile Dressing for Pressure Relief and Hygiene:

LIGASANO® white unsterile is suitable for single use only. The product must not be used after the expiration date and is not suitable for reconditioning. LIGASANO® white unsterile must not continuously remain on the body. Recommendation: dressing change after 1-3 days, depending upon the indication. Any applications in contact, connection or combination with any additional preparations, drugs, solutions, ointments etc., are not guaranteed by us. Applications in combination with additional mechanical, electrical or electronical apparatuses or aids are not guaranteed by us.

► **Prevention of decubitus ulcers through active positioning of the patient:** the rapid drop in pressure tension from LIGASANO® white enables almost even distribution of the pressure from the surface padding and thus avoids the dreaded pressure points. LIGASANO® white can be simply cut to fit. This is how aids such as pulleys, wedges and heel shoes etc. are created for changing patients' positions, having been adapted by imaginative nursing staff. The mechanical stimulus of LIGASANO® white stimulates the peripheral blood supply where there is contact with the skin, and the patient's residual pressure increases. Surplus sweat is absorbed, leaving normal skin moisture, and care of the skin is also associated with the mechanical stimulus. If LIGASANO® white is the appropriate thickness it does not form any folds, thus actively preventing a further risk of decubitus. **Patients are always placed directly onto LIGASANO® white with no layer in between so that there is skin contact.**

DE

GB

FR

IT

ES

RO

SE

NL

RUS

PL

CZ

LT

GR

HU

TR

► **Prevention against fungal infections and/or inflammations in skin-skin-contact:** a wet and warm low-oxygen environment around the skin /skin-contact encourages the occurrence of fungal infections and inflammations of the skin. An intermediate layer of 1-2 cm of LIGASANO® white is extremely effective against this, being both preventive and therapeutic. If fungal infections do occur, an additional remedy can be applied thinly to the affected sections of skin, if necessary. Additional measures are usually not necessary in case of local inflammation. **Typical application:** Between toes and fingers or in folds of the belly, under breasts, and for armpits and contractures.

Sterilisation on-site and subsequent application of the thus sterilised products for wound treatment:

The non-sterile product has to be sterilised at 134 °C for five minutes with moist heat according to a validated procedure. These standardised parameters meet the requirements for adequate product safety and economy, and there are no measurable effects on the properties of LIGASANO® white.

When sterilising unsterile LIGASANO® products on site (e.g. in the clinic, practice, etc.), the legal regulations and risks must be observed. In particular, we point out that if the contamination is too high (> 200 CFU), the sterilised product could be contaminated with an excessive, risk-increasing amount of pathogenic germs.

In case of sterilisation or resterilisation you act on your own responsibility. Please also note that the indication, side effects, contraindication, performance characteristics and use change when you apply the on-site sterilised product for wound treatment.

Intended Use for On-site Sterilised Products:

Promotion of wound healing in sterile application by cleaning, covering or filling wounds.

Indications for On-site Sterilised Products:

Acute and chronic wounds (e.g. decubital ulcers, ulcer cruris, diabetic foot ulcers). Post-surgical wounds and post-surgical wound healing impairments. Thermal wounds, such as burns, scalds, frostbite.

Adverse Effects, Side Effects, Residual Risks for On-site Sterilised Products:

Wound cleaning: mechanical debridement may cause pain, even with LIGASANO® white sterile.

Wound treatment: because of its cell structure a vascularisation (e.g. of granulation tissue) into the material is possible. You can prevent this with a timely dressing change, after 1-3 days using LIGASANO® white sterile. On passive wounds or parts of the body that receive poor blood supply, pain may occur after application. This occurs when the blood circulation is stimulated by the mechanical stimulus of the material to such an extent, that previously diminished or reduced pain returns. Initially the pain can be intense, however it is usually regulated after a few hours to days. If this effect is not desired, you should forgo its application. Too much pressure due to the material's prior tension, external influence, or suction may lead to compression of vessels and therefore to pressure ulcers.

It may occur that LIGASANO® white sterile sticks to the wound, mostly due to too much air reaching the wound. This means that the LIGASANO® white sterile dressing is not thick enough (it should be 2 cm from skin level). LIGASANO® white can be moistened on the side facing the wound (e.g. with Ringer's solution or wound irrigation solution) to prevent it from sticking.

If hypergranulation occurs, please treat it appropriately and switch to other dressing material groups, if necessary.

LIGASANO® white is a very simple product and its effects are due to its structure. It does not contain nor emit any active compounds. With the correct application of LIGASANO® white there are no known negative side effects, nor incompatibility or interactions with medicine. Typical side effects may include reddening of the skin upon contact and initial "itching". More intense side effects may include return of sensitivity due to stimulation of local blood circulation. Its effects are only physical and end almost immediately after discontinuation of application. As with all wound dressings, minor skin reactions such as maceration, erythema, secondary infected dermatitis, erysipelas as well as hypersensitivity reactions and pain at dressing change may occur.

Please report any incidents that have occurred in connection with our dressing material LIGASANO® white sterile to us or to the competent authority of the member state in which the incident occurred.

Contraindications for On-site Sterilised Products:

Tumour wounds and untreated osteomyelitis. Not suitable for contact with organs, contact to exposed blood vessels with the danger of damage, or contact with nerves.

Patient Target Group for On-site Sterilised Products:

Children (> 2 - 16 years of age) and adults (> 16 years of age) with acute, chronic, post-surgical wounds/wound healing disturbances and thermal wounds.

User/Practitioner of the Products Sterilised On-site:

Application must always be done by healthcare professionals. The decision about the particular kind of application of LIGASANO® white sterile is determined by healthcare professionals and is adapted and appropriated in every single case.

Capability Characteristics and Clinical Benefit of the Products Sterilised On-site:

With the sterile wound dressing LIGASANO® white you can **clean, cover or pack wounds**. The product can be used both for **mechanical debridement** and for **wound treatment in all phases of wound healing**. The expected clinical benefit presents itself in a very effective wound cleaning, wound bed conditioning and in the healing of the wounds.

Application of On-Site Sterilised LIGASANO® white for Wound Treatment:

LIGASANO® white sterile is suitable for single use only. Sterility is ensured only by intact packaging. The product must not be used after the expiration date and is not suitable for reconditioning or resterilisation due to its potential for contamination.

LIGASANO® white sterile must not continuously remain in the wound or on the body. The application may be repeated over a period up to 12 months, always with a new dressing. Dressing change after 1-3 days, depending upon the indication. Please consider that LIGASANO® white sterile must not have direct contact to organs. It may have direct contact to mucosa.

Any applications in contact, connection or combination with any additional preparations, drugs, solutions, ointments etc., are not guaranteed by us. Applications in combination with additional mechanical, electrical or electronical apparatuses or aids are not guaranteed by us.

► Wound Cleaning / Mechanical Debridement:

LIGASANO® white sterile is suited for mechanical cleaning of nearly every kind of wound (see contraindications) and during all wound phases. It is highly effective in removing soft coatings and biofilm. The wound and the edge of the wound can be treated with LIGASANO® white sterile. For this treatment LIGASANO® white sterile remains in the wound only for a short time to clean it. You can treat the wound with either the twisting or the wiping technique. Please ensure that the wiping direction must always be from inside to outside. Important: please note that LIGASANO® used for wound cleaning must not remain in the wound. If the proceeding wound dressing is to be made with LIGASANO®, a new wound pad must be applied.

► Wound Treatment:

LIGASANO® white sterile is usually applied directly onto or into the wound without any further steps or products (e.g. ointments or irrigation). Additional measures may negatively affect the desired result. Apply a 1-2cm thick LIGASANO® white sterile to skin-surface wounds, the edges overlapping the wound by 1-2cm. Deeper wounds, pocket-shaped wounds and hollow wounds should be **completely packed with sterile LIGASANO® white**. There must be reliable contact completely covering the wound. To guarantee this, LIGASANO® white sterile can be compressed.

Padding thickness: measuring from the wound base, LIGASANO® white sterile should be at least 2cm thick. The thickness can also be formed using several layers, including the outer covering of the wound. If the wound dressing composition consists of several sheets or pieces, please count the number of used dressing pieces and note it in the wound documentation, to ensure the removal of all applied pieces.

Dressing change frequency: one time every 1-3 days, based on necessity, and as soon as the first sign of secretion becomes visible on the outside of LIGASANO® white sterile. For wounds that heavily suppurate during the cleaning phase, the dressing should also be changed several times daily.

Padding technique: each case should be evaluated and adapted accordingly. Ensure that during dressing change, the wound is left uncovered for the shortest possible amount of time in order to avoid cooling down and to keep risk of re-infection low.

Typical wound cleaning sequence: with chronic wounds the wound-cleaning phase requires approx. 10-20 dressing changes. The intended result is the smallest possible "clean" wound, with low susceptibility to re-infection.

Storage, Transport, Shelf Life:


Store LIGASANO® white unsterile in a dry area that is protected against UV Radiation. The original packaging offers protection against humidity and UV radiation due to the filtering plastic film, and therefore storage is only permitted in the original packaging. The product may display a yellowish discolouration, especially if improperly stored, however this does not affect the quality. Special transport conditions (e.g. cooling) are not necessary. The product is stable for 36 months after the date of manufacture. The date of manufacture and expiration can be found on the product label of the original packaging.

Disposal / Environmental Relevance:

In addition to its accredited quality (QM system for medical devices), LIGAMED® medical Produkte GmbH also documents environmental awareness by validation according to EC regulation and its participation in „Umweltpakt Bayern“ (Bavarian Environmental Understanding). We will send you our current certificates upon request. All LIGASANO® products are single-substance-materials made of polyurethane,

consisting of approx. 97% air, and can be disposed of with ordinary household waste. Disposal of the packaging is internationally regulated in different ways. In Germany it is regulated by a licensed company.

Symbols Used for Labelling:

	Item number		Manufacturer
	Lot number		Production date
	Use by: year / month		Distributor
	Medical device		Unsterile
	Unique device identifier		Disposable product
	Consult instructions for use		Protect from sun and UV radiation
	Do not use if the packaging is damaged		Protect from moisture

Date of Revision:

2023 / 07

Information about the Manufacturer:

LIGAMED® medical Produkte GmbH • Pfannenstielstr. 2 • 90556 Cadolzburg/Germany • Tel. +49 9103 2046 • E-Mail info@ligamed.de

An electronic version of this instruction for use is available on our website: www.ligasano.com/downloads

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A.Pf.E.L.®,  are Trade Marks of LIGAMED® medical Produkte GmbH, Germany
Copyright by LIGAMED® medical Produkte GmbH
LIGAMED® medical Produkte GmbH • Pfannenstielstr. 2 • D - 90556 Cadolzburg