REVMIL Collagenverband

NL

SAMENSTELLING

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BiologQ, Apeldoorn (NL) – Bfatory, Rhenen (NL).

REVMIL

Wound medication based on collagen and pure honey

EN

COMPOSITION

Each sterile pad is made up of type I equine collagen (60%) and pure honey (40%).

DESCRIPTION

This medical device is a golden brown, sterile spongy pad containing type I heterologous collagen extracted from equine tendons, with the addition of pure honey. The pads are easily adaptable to different body areas; upon application onto wounds, these collagen- and honey-based products create a protective layer that fosters natural healing processes and enhances antibacterial protection. Collagen acts as a scaffold for fibroblasts, favors the restoration of optimal connective tissue conditions and provides the necessary support to natural wound regeneration. Honey helps maintaining a moist wound environment and contributes to protecting lesions against bacterial infections.

INDICATIONS

This medical device is indicated for the following: treatment of cutaneous lesions of different etiology, including acute and chronic delayed-healing ulcers, venous ulcers, diabetic ulcers, pressure sores, infected wounds, cancer wounds, as well as first and second degree burn wounds; in case of vascular reconstructive surgery, vascular surgery, carotid surgery, abdominal and gynecological surgery, in orthopedics, traumatology and dentistry as a support to physiological wound healing processes; local hemostasis, to control capillary bleeding during first aid.

CONTRAINDICATIONS

Do not use the device on the conjunctiva. Do not use it on patients with a family history of known anaphylactic reactions, allergies or individual hypersensitivity to collagen or honey.

METHODS OF USE AND POSOLOGY

Use the product immediately after opening the pad's primary packaging. When the blister has been opened, if only part of the pad is used, do not reuse any remaining parts of the pad for subsequent applications. Before using the product on infected wounds, perform a systemic antibiotic treatment and/or a topical antimicrobial treatment.

Pre-treat the wound bed and clean it by repeatedly dabbing it with wet gauze soaked in saline solution. After debriding the wound, apply the pad on it. If exudate is insufficient, moistenize the pad with saline solution. To cover the whole lesion and ensure firm, full contact between the whole area to be treated and the pad, several pads can be placed one after the other or overlapped. If needed, fasten the pad to the lesion with sterile gauze or non-adherent dressing, then apply an elastic compression bandage.

When the device is placed upon the exudate, it perfectly adapts to the wound; then it is absorbed by the body and it dissolves after some time.

It is recommended to apply the dressing every 3-4 days or more frequently, depending on the wound exudate. The wound should be monitored every 2-3 days and the treatment should be continued until complete healing. If the pad has not been absorbed when a new application or a wound follow-up visit is due, do not remove the pad from the wound bed.

If the device was used as a hemostatic agent, and there is no certainty as to whether application was executed under aseptic conditions, remove the device when hemostasis is achieved.

WARNINGS AND PRECAUTIONS

Despite its hemostatic properties, the device is not to be used instead of a tourniquet or direct compression procedures in case of severe bleeding. Make sure that the package and the blister are intact before use. Sterility and quality of the device can be guaranteed only for blisters/packages sealed at the origin.

Always the product or using it in any way that is not in accordance with the conditions of use described in this document may affect its sterility and performance, which will not be guaranteed any longer. The device is sterile. Do not use it if the package is damaged.

The device is for single use; do not reuse any remaining parts of a pad after its first opening. Do not use the product after its expiry date. The expiry date refers to correctly stored products with intact packaging. There are no known contraindications to the use of the product during pregnancy or breast-feeding; however, in the absence of specific data, it is recommended to use it under direct medical supervision.

The product is not toxic if ingested; it does not form a bolus; it does not swell with moisture, therefore it does not pose any choking risk. In case of unexpected side effects, seek medical advice.

The device will not affect the ability to drive vehicles or operate machinery. *KEEP THE DEVICE OUT OF THE SIGHT AND REACH OF CHILDREN.*

INTERACTION

There are no known interactions with other substances. The dressings' activity is not affected by anticoagulant therapies.

UNDESIRABLE SIDE EFFECTS

The device is biocompatible and bioresorbable. No cases of sensitization or undesirable side effects have been observed or reported for its intended use. In case of side effects, seek medical advice.

Any side effects should be reported to the distributor and /or manufacturer.

TRANSPORTATION AND STORAGE CONDITIONS

Revmil must be transported and stored at a temperature not exceeding 25 °C; it must be kept in a cool, dry place away from heat sources, avoiding exposure to direct sunlight. Do not keep the product after the expiry date indicated on the package.

CONTENTS OF THE PACKAGE

Cardboard box containing 3 5x5 cm sterile pads in individual blisters.

MANUFACTURER

Eurosearch s.r.l. - Via Larga, 15 – 20122 Milan (I) - Ph. +39 02 8055660.

DISTRIBUTOR

Oswell Penda Pharmaceuticals, Oswestry Shropshire (UK).

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