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### Product

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## Patient Information Leaflet Revamil®

### The product

Revamil® single dose, hydrophilic wound gel  
 2 g honey for the treatment of wounds

### Product description

Revamil® is a wound gel that consists of pure honey, which has been produced under controlled conditions. The carefully controlled production process guarantees that Revamil® does not contain any traces of pesticides and that the same quality can be reproduced consistently. The honey used to make Revamil® contains a high proportion of honey enzymes such as glucose oxidase. The highly concentrated sugars cause some of the wound exudate to mix with Revamil® gel. The resulting, covering layer helps to strengthen the antibacterial protection provided by the dressing. The product has a cleansing function as it creates a moist wound environment which stimulates debridement in the wound bed.

### Instructions for use

1. Before applying Revamil® gel first clean the wound.
2. Remove the closing tip from the syringe.
3. By slowly pressing the plunger of the syringe Revamil® wound gel is released from the syringe.
4. Spread a thin layer of Revamil® on the secondary wound dressing.
5. Revamil® can also be applied directly to the wound and then covered with a secondary dressing.
6. Attach the secondary dressing in the normal way.
7. Depending on the nature of the wound, replace the dressing every day or every two days.

### Indications

Revamil® wound gel can be used to treat chronic sores (bed sores and different kinds of ulcers), infected wounds, acute wounds, surgical oncological wounds, radiation induced oncological wounds and minor (first- and second-degree) burns.

### Contra indications

- Do not use the product if the patient is known to be allergic or hypersensitive to honey. An allergic reaction or hypersensitivity can cause redness, irritation and possibly minor pain in the skin or wound.

### Side effects

- During the first minutes after application, the product can cause a light stinging sensation. The stinging sensation is probably caused by the low pH and high osmolarity of the honey and has no effect on the performance of the product.
- In the event of unexpected side effects after use, contact your doctor or the manufacturer. Side effects can cause inconvenience to the patient.

### Safety precautions

- Do not use if honey is leaking from the inner packaging. Leakage can lead to contamination of the product with bacteria and to wound infection. Because of the antibacterial protection that honey provides, the chance of wound infection is considered very low.
- Do not use if the syringe is damaged as this can lead to contamination of the product with bacteria and to wound infection. Because of the antibacterial protection that honey provides, the chance of wound infection is considered very low.
- Do not allow the syringe to come into direct contact with the surface of the wound. Direct contact between the wound/skin and the uncapped syringe could lead to contamination with bacteria and wound infection. Because of the antibacterial protection that honey provides, the chance of wound infection is considered very low.
- Before applying Revamil®, dampen the wound with 0.9% NaCl. Application of Revamil® Gel on a dry wound could lead to lower performance of the product and slower wound healing.
- Do not apply this product for more than 30 days continuously. Treatment for more than 30 days continuously could lead to inconvenience to the patient, such as irritation and/or minor pain.
- Revamil® Single Dose provides only one dose per unit. Under-dosage of the product could lead to lower performance of the product and slower wound healing.
- Do not treat different patients with the same syringe as this can lead to cross-contamination with bacteria and to wound infection. Because of the antibacterial protection that honey provides, the chance of wound infection is considered very low.
- Do not gamma irradiate a Revamil® Single Dose syringe again. Re-gamma irradiation of an already used syringe could lead to contamination of the product and to wound infection.
- Do not reuse an already used syringe. Reuse of a syringe can lead to loss of function of the honey, contamination of the product and/or wound infection.
- Do not use the product after the expiration date as this can lead to lower performance of the product and to inconvenience for the patient.

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Bfactory Health Products B.V.
DK: Axel Madsen A/S, Vedbaek
DE: ROGG Verbandstoffe GmbH & Co. KG, Fahrenzhausen
FI: Latua Group Oy, Helsinki