Distributed by

**puras**

031 996 85 85  
puras.ch

**OPERATOR'S MANUAL**  
**Pulse Oximeter ri-fox N**

ENGLISH

## ***General Description***

Oxygen Saturation is a very important parameter for the Respiratory Circulation System. Many respiratory diseases can result in oxygen saturation being lowered in human blood. Additionally, the following factors can reduce oxygen saturation: organ dysfunction caused by anesthesia, intensive postoperative trauma, injuries caused by some medical examinations. That situation might result in light-headedness, asthenia, and vomiting. Therefore, it is very important to know the oxygen saturation of a patient so that doctors can find problems in a timely manner.

The fingertip pulse oximeter ri-fox N features small size, low power consumption, convenient operation and portability. It is only necessary for a patient to put one of his fingers into the fingertip photoelectric sensor for diagnosis, and a display screen will show oxygen saturation. It has been proven in clinical experiments that it also features high precision and repeatability.

## ***Product Operation Scope***

ri-fox N is a portable non-invasive, spot-check, oxygen saturation of arterial hemoglobin (SpO<sub>2</sub>) and pulse rate of adult and pediatric patients at home and hospital (including clinical use in internal medicine/surgery, anesthesia, etc). The device is not intended for continuous monitoring.

The ri-fox N requires no routine calibration or maintenance other than replacement of batteries.

## ***Technical Specifications***

**1. Display:** LED

**2. SpO<sub>2</sub>:**

Measurement range: 70~100%

Accuracy: 70% ~ 100%,  $\pm 2\%$ ; <70% no definition

**3. Pulse Rate:**

Measure range: 30~235 bpm

Accuracy: 30~100bpm  $\pm 2$ bpm; 101~235bpm,  $\pm 2\%$

Pulse Intensity: Bar graph Indicator

**4. Power Requirements:**

Two AAA alkaline Batteries

Power consumption: Less than 30mA

Low power indication:

Battery Life: Two AAA 1.5V, 800mAh alkaline batteries could be continuously operated as long as 30 hours.

**5. Dimension:**

Length: 66mm

Width: 39mm

Height: 32mm

Weight: 34g (without batteries)

**6. Environment Requirements:**

Operation Temperature: 5°C~40°C

Storage Temperature: -20~55°C

Ambient Humidity: ≤80%, no condensation in operation;  
≤93%, no condensation in storage

**7. Measurement Performance in Low Perfusion Condition:** The pulse wave is available without failure when the simulation pulse wave amplitude is at 0.6% using the test equipment (Biotek Index Pulse Oximeter Tester).

**8. Interference Resistance Capacity against Ambient Light:** Device works normally when mixed noise produced by Biotek Index Pulse Oximeter Tester.

## ***Precautions for use***

1. Do not use the pulse oximeter in an MRI or CT environment
2. Do not use the pulse oximeter in situations where alarms are required. The device has no alarms.
3. Explosion hazard: Do not use the pulse oximeter in an explosive atmosphere.
4. The pulse oximeter is intended only as an adjunct in patient assessment. It must be used in conjunction with other methods of assessing clinical signs and symptoms.
5. Check the pulse oximeter sensor application site frequently to determine the positioning of the sensor and circulation and skin sensitivity of the patient.
6. Before use, carefully read the manual.
7. ri-fox N has no SpO<sub>2</sub> alarms; it is not for continuous monitoring, as indicated by the symbol.
8. Prolonged use or the patient's condition may require changing the sensor site periodically. Change sensor site and check skin integrity, circulatory status, and correct alignment at least every 4 hours.

9. Inaccurate measurements may be caused by autoclaving, ethylene oxide sterilizing, or immersing the sensors in liquid may cause inaccurate readings.
10. SpO<sub>2</sub> measurements may be adversely affected in the presence of high ambient light. Shield the sensor area if necessary.
11. Excessive patient movement may cause inaccurate readings.
12. Placement of a sensor on an extremity with a blood pressure cuff, arterial catheter, or intravascular line may cause inaccurate readings.
13. Fingernail polish or false fingernails may cause inaccurate SpO<sub>2</sub> readings.

#### **Contraindication**

14. Significant levels of dysfunctional hemoglobins (such as carbonxy-hemoglobin or methemoglobin) may cause inaccurate readings
15. Intravascular dyes such as indocyanine green or methylene blue may cause inaccurate readings
16. Venous pulsations may cause inaccurate readings.
17. Hypotension, severe vasoconstriction, severe anemia or hypothermia may cause inaccurate readings.
18. Cardiac arrest or shock may cause inaccurate readings.

#### **No side effects known**

***Follow local ordinances and recycling instructions regarding disposal or recycling of the device and device components, including batteries.***


## ***Product specifications***

1. Easy to use.
2. Small in volume, light in weight.
3. Low power consumption.
4. Low voltage warning will be indicated in visual window when battery voltage is low.
5. The product will automatically be powered off when no signal is registered for longer than 8 seconds.


## ***Operation Instructions***

1. Install two AAA batteries into the battery cassette, and then close its cover.
2. Open the clamp as shown in the picture below.
3. Insert one finger completely into the opening of the pulse oximeter.
4. Press the power switch once on the front panel.
5. The patient's finger should remain still while the pulse oximeter is working.
6. Read corresponding data from display screen.

***When you insert a finger into the Oximeter, your nail surface must be upward.***

## ***Brief Description of Front Panel***

The Pulse Bar graph displays the strength of the pulse rate signal. The height of the bar graph shows the patient's pulse strength.

## ***Battery Installation***

1. Put the two AAA batteries into battery cassette in correct polarities.
2. Press the button down along the arrow and push the battery cover shown as below:

**Press the button down!**

## ***Maintenance and Storage***

1. Replace batteries in time when low voltage lamp is lighted.
2. Clean surface of the ri-fox N before it is used in diagnosis for patients.
3. Remove batteries inside the battery cassette if the oximeter will not be used for a long time.
4. It is best to preserve the product in a place where ambient temperature is  $-20\sim 55^{\circ}\text{C}$  and humidity is  $<93\%$ , no condensation.
5. It is recommended that the product should be kept in a dry environment anytime. A wet ambient might affect its lifetime and even might damage the product.
6. Please follow the law of the local government to deal with used batteries.

## ***Calibrating the Oximeter***

1. The functional tester cannot be used to assess the accuracy of the oximeter.
2. The test method used to establish the  $\text{SpO}_2$  accuracy is clinical testing. The oximeters used to measure the arterial haemoglobin oxygen saturation levels are compared to the levels determined from arterial blood sampling with a CO-oximeter.
3. Index 2 that is made by Biotek company is a function tester. Set Tech to 1, R curve to 2, and then users can use this particular calibration curve to measure the oximeter.

## ***Declaration***

EMC of this product complies with IEC60601-1-2 standard. Portable and mobile high frequency appliances can influence medical electronic devices. Medical electronic devices need special safety requirements concerning EMC. **You may find a detailed description of the manufacturers EMC declaration at the end of this user manual.**

The materials which the user can come into contact have no toxicity and no action on tissues, comply with ISO10993-1, ISO10993-5 and ISO10993-10.

### **Guidance and manufacture's declaration – electromagnetic emissions-for all EQUIPMENT and SYSTEMS**

Guidance and manufacture's declaration – electromagnetic emission		
The ri-fox N is intended for use in the electromagnetic environment specified below. The customer or the user of the ri-fox N should assure that it is used in such and environment.		
<b>Emission test</b>	<b>Compliance</b>	<b>Electromagnetic environment – guidance</b>
RF emissions CISPR 11	Group 1	The ri-fox N uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emission CISPR 11	Class B	The ri-fox N is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.

## Possible Problems and resolutions

Problems	Possible reason	Solution
SpO <sub>2</sub> or PR can not be shown normally	<ol style="list-style-type: none"> <li>1.Finger is not plugged correctly</li> <li>2.Patient's Oxyhemoglobin value is too low to be measured.</li> </ol>	<ol style="list-style-type: none"> <li>1. Retry by plugging the finger.</li> <li>2. Try some more times, If you can make sure about no problem existing in the product. Please go to a hospital for exact diagnosis.</li> </ol>
SpO <sub>2</sub> or PR is shown unstably	<ol style="list-style-type: none"> <li>1.Finger might not be plugged deeply enough.</li> <li>2.Finger is trembling or patient's body is in movement status.</li> </ol>	<ol style="list-style-type: none"> <li>1.Retry by plugging the finger.</li> <li>2.Try not to move.</li> </ol>
The Oximeter can not be powered on	<ol style="list-style-type: none"> <li>1.Power of batteries might be inadequate or not be there at all.</li> <li>2.Batteries might be installed incorrectly.</li> <li>3.The Oximeter might be damaged.</li> </ol>	<ol style="list-style-type: none"> <li>1.Please replace batteries</li> <li>2.Please reinstall the batteries</li> <li>3.Please contact with local customer service centre.</li> </ol>
Indication lamps are suddenly off	<ol style="list-style-type: none"> <li>1.The product is automatically powered off when no signal is detected longer than 8 seconds</li> <li>2. Lower power</li> </ol>	<ol style="list-style-type: none"> <li>1.Normal</li> <li>2.Replace the batteries</li> </ol>
"Error3" or "Error4" is displayed on screen.	<ol style="list-style-type: none"> <li>1.Low power</li> <li>2.Mechanical malfunction for receive-emission tube</li> <li>3. Amp circuit malfunction.</li> </ol>	<ol style="list-style-type: none"> <li>1.Change batteries.</li> <li>2.Please contact the local customer service center.</li> <li>3.Please contact the local customer service center.</li> </ol>

## Symbol Definitions

Symbol	Definition
	Type BF applied part
	Attention, consult accompanying documents.
%SpO <sub>2</sub>	Oxygen saturation
♥ PR BPM	Pulse rate (BPM)
	Low power indication
	NOT for continuous monitoring
SN	Serial No.

**Note:**

*The illustrations used in this manual may differ slightly from the appearance of the actual product.*

## ***Disinfection***

Please use medical alcohol to clean the plastic material touching the finger inside of the pulse oximeter, and clean the test sensor using alcohol before and after operation.

## ***Warranty***

This product has been manufactured under the strictest quality standards and has undergone a thorough final quality check before leaving our factory. We are therefore pleased to be able to provide a warranty of 30 months from the date of purchase on all defects, which can verifiably be shown to be due to material or manufacturing faults.

A warranty claim does not apply in the case of improper handling. All defective parts of the product will be replaced or repaired free of charge within the warranty period. This does not apply to wearing parts.

Please remember that all warranty claims have to be made during the warranty period. We will, of course, be pleased to carry out checks or repairs after expiry of the warranty period at a charge. You are also welcome to request a provisional cost estimate from us free of charge. In case of a warranty claim or repair, please return the RIESTER product to your dealer.

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**Guidance and manufacturer's declaration – electromagnetic emission –  
for all EQUIPMENT AND SYSTEMS**

Row

1	Guidance and manufacturer's declaration – electromagnetic emission		
2	The MD300C61 and MD300C63 Pulse Oximeter are intended for use in the electromagnetic environment specified below. The customer or the user of MD300C61 and MD300C63 Pulse Oximeter should assure that it is used in such an environment.		
3	Emissions test	Compliance	Electromagnetic environment - guidance
4	RF emissions CISPR 11	Group 1	The MD300C61 and MD300C63 Pulse Oximeter use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
5	RF emissions CISPR 11	Class B	
6	Harmonic emissions IEC 61000-3-2	N/A	
7	Voltage fluctuations / flicker emissions IEC 61000-3-3	N/A	

**Guidance and manufacturer's declaration – electromagnetic immunity –  
for all EQUIPMENT and SYSTEMS**

<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The MD300C61 and MD300C63 Pulse Oximeter are intended for use in the electromagnetic environment specified below. The customer or the user of the MD300C61 and MD300C63 Pulse Oximeter should assure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 test level</b>	<b>Compliance level</b>	<b>Electromagnetic environment - guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrostatic transient / burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	± 1 kV differential mode ± 2 kV common mode	N/A	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	< 5 % $U_T$ (-95 % dip in $U_T$ ) for 0.5 cycle  40 % $U_T$ (60 % dip in $U_T$ ) for 5 cycles  70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles  < 5 % $U_T$ (-95 % dip in $U_T$ ) for 5 sec	N/A	Mains power quality should be that of a typical commercial or hospital environment. If the user of the MD300C61 and MD300C63 Pulse Oximeter requires continued operation during power mains interruptions, it is recommended that the MD300C61 and MD300C63 Pulse Oximeter be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: $U_T$ is the a. c. mains voltage prior to application of the test level.			

**Guidance and manufacturer's declaration – electromagnetic immunity –  
for EQUIPMENT and SYSTEM that are not LIFE-SUPPORTING**

Guidance and manufacturer's declaration – electromagnetic immunity			
The MD300C61 and MD300C63 Pulse Oximeter are intended for use in the electromagnetic environment specified below. The customer or the user of the MD300C61 and MD300C63 Pulse Oximeter should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	N/A	<p>Portable and mobile RF communications equipment should be used no closer to any part of the MD300C61 and MD300C63 Pulse Oximeter, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter:</p> <p><b>Recommended separation distance</b></p> $d = \left[ \frac{3.5}{P} \right] \sqrt{P}$ $d = \left[ \frac{3.5}{E} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[ \frac{7}{E} \right] \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in metres (m).<sup>1</sup></p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>2</sup> should be less than the compliance level in each frequency range.<sup>3</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.			
NOTE 2: These guidelines may not apply in all situations. Electromagnetic is affected by absorption and reflection from structures, objects and people.			
<p><sup>1</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio-broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the MD300C61 and MD300C63 Pulse Oximeter is used exceeds the applicable RF compliance level above, the MD300C61 and MD300C63 Pulse Oximeter should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the MD300C61 and MD300C63 Pulse Oximeter.</p> <p><sup>2</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.</p>			

**Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM - for EQUIPMENT and SYSTEMS that are not LIFE-SUPPORTING**

Recommended separation distances between portable and mobile RF communications equipment and the MD300C61 and MD300C63 Pulse Oximeter			
The MD300C61 and MD300C63 Pulse Oximeter are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the MD300C61 and MD300C63 Pulse Oximeter can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the MD300C61 and MD300C63 Pulse Oximeter as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[ \frac{3.5}{F_1} \right] \sqrt{P}$	80 MHz to 800 MHz $d = \left[ \frac{3.5}{E_1} \right] \sqrt{P}$	800 MHz to 2.5 GHz $d = \left[ \frac{7}{E_1} \right] \sqrt{P}$
0.01	/	0.12	0.23
0.1	/	0.38	0.73
1	/	1.2	2.3
10	/	3.8	7.3
100	/	12	23
For transmitters rated at a maximum output power not listed above the recommended separation distance <i>d</i> in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where <i>P</i> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.			
NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			