Instructions for Use

GB	Great Britain		_____	P1
DE	German	Deutsch	_____	P2
ES	Spanish	Español	_____	P3
FR	French	Français	_____	P4
IT	Italian	Italiano	_____	P5
LV	Latvian	Latviski	_____	P6
NL	Dutch	Nederlands	_____	P7

Effective

Effective

<Instructions for Use>

Please read these instructions carefully before use.

Indications

Winged Needle sets are intended for use as part of a device for introducing solutions or collecting blood from the patient's vascular system and are designed to connect to infusion sets, transfusion sets, blood collection sets and others.

Contraindications

The use of winged needle sets is not in itself a therapeutic procedure and there are no particular contraindications when they are used for the intended purposes.

Instructions for Use

1. Remove winged needle set from packaging.
2. Check to make sure that safety mechanism is securely locked into the stopper. (Figure 1)
3. Remove connector cap from luer connector, connect the luer connector to the infusion line and then prime the line (Figure 2), or connect the luer connector to the blood collection set.
4. Remove needle cap by firmly holding wings and being cautious with the exposed needle tip.
5. Insert needle into the disinfected puncture site.
6. Make sure that the needle tip is securely in the blood vessel by confirming blood flashback visibly in the safety mechanism and immobilize the wings by taping, if required.
7. For infusion, adjust the optional infusion set clamp, make sure that the medication is flowing, then begin the infusion.
8. For blood collection, by attaching a luer adaptor, you can collect specimen into vacuum blood collection tubes.

Activation of Safety Mechanism

The following procedures are used to activate the accidental needlestick prevention mechanism. We recommend using a pre-removal activation.

Pre-Removal Activation

1. After completion of the infusion or collection, close the optional pad and apply digital pressure to the puncture site using sterilized gauze pad and, using the same hand, hold down the safety mechanism and/or the wings. (Figure 3)
2. With the other hand, release the lock by pressing in both sides of the stopper. (Figure 3)
3. NOTE : Do not apply undue pressure to the puncture site, slide the stopper backwards until it becomes completely locked with an audible or tactile click. (Figure 3)
4. Verify that the stopper is securely immobilized against the safety mechanism. Therefore, the needle is completely retracted and protected by the safety mechanism. (Figure 4)
5. Maintain pressure on the puncture site according to facility protocol.
6. If used, remove the tape from the winged needle set and promptly dispose of the set. (See Cautions 2)

Warnings






1. The device will perform correctly if the instructions are followed appropriately. Please refrain from carelessly releasing the lock, forcibly pulling the wings or doing any such actions, since rough handling of the device may damage it.
2. Due to the risk of coming into contact with blood, please wear gloves and eye protector.
3. Do not forcibly release or destroy the accidental needlestick prevention mechanism after it has been activated.

Cautions

1. When immobilizing the wings, be careful NOT TO COVER THE STOPPER WITH TAPE since taping the stopper prevents activating the safety mechanism. (Figure 5)
2. Due to the risk of blood dripping from its exposed tip, promptly dispose of the product pointing the needle tip upwards in an approved disposal container after activating the safety mechanism.
3. When the safety mechanism is not being activated, keep hands away from the needle tip at all times.
4. There is a possibility of blood, medicine, etc., leaking out of the product when removing the luer connector. Be careful.
5. Immediately use the product after breaking package seal.
6. In the event of damaged packaging, do not use the product.
7. Store product in a cool, dry environment away from direct sunlight, extreme temperature and humidity.
8. While using and during post-use disposal, pay close attention to prevent accidental needlestick.
9. This product is for single use only. DO NOT REUSE.
Reuse or reprocessing of a single use device may lead to contamination and compromised device function or structural integrity
10. The use of an alcohol or solvent-based solution may adversely affect the device tightness (connector integrity) leading to leaks. If required, use only aqueous-based solution such as povidone-iodine solution, (i.e., Betadine®)
11. Be careful not to excessively pull, push into, or put load that may bent the joint between tube and connector.
12. Be careful not to damage the tube with sharps such as scissors, or by picking up with forceps, etc., since it may damage the tube or leakage.
13. Be careful not to put excessive load on needle, since it may bend or break the needle, or damage the hub.

GUARANTEE

1. Our winged needle sets are manufactured under strict quality control and the quality is assured. We shall not be responsible for injury caused to a patient or any person or damage to any object that is attributed to transport, storage, and usage in your institution.
2. If a patient or any person is injured or any object is damaged when using our winged needle sets, we shall not be responsible for the injury or damage unless we are clearly identified as being at fault.
3. If a patient or any person is injured or any object is damaged when reusing our winged needle sets, we shall not be responsible for the injury or damage of any nature.
4. We shall not be responsible for any injury or damage caused if our winged needle sets are used after the expiry date mentioned on the packages.

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




Effective

Effective

	Do not re-use ¹		Manufacturer ⁷
	Consult instructions for use ²		Date of manufacture ⁸
	Sterilized using ethylene oxide ³		Keep dry ⁹
	Do not use if package is damaged ⁴		Keep away from sunlight ¹⁰
	Use-by date ⁵		Not made with natural rubber latex ¹¹
	Batch code ⁶		NON-Bis(2-ethylhexyl) phthalate ¹²

GB Symbols used for labeling

- 1 Do not re-use
- 2 Consult instructions for use
- 3 Sterilized using ethylene oxide
- 4 Do not use if package is damaged
- 5 Use-by date
- 6 Batch code
- 7 Manufacturer
- 8 Date of manufacture
- 9 Keep dry
- 10 Keep away from sunlight
- 11 Not made with natural rubber latex
- 12 NON-Bis(2-ethylhexyl) phthalate

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Title: B. Braun Venofix Safety IFU.pdf Initiator: Stefan ? Klaas

This document is signed electronically in compliance with the B. Braun electronic signature policies and procedures by following persons:

UserName: Klaas, Stefan (klaastde)
Title: Product Manager IV-Catheters & Scalp Vein Sets
Date: Monday, 31 October 2022, 21:09 W. Europe Daylight Time
Meaning: Document signed as Author

UserName: Padberg, Andreas (padbande)
Title: HC CoE Infusion and Pain Therapy - Manager Regulatory Affairs
Date: Wednesday, 02 November 2022, 08:15 W. Europe Daylight Time
Meaning: Approve Document

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Nipro Corporation
B.Braun Melsungen AG