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CE CE 2460 - DNV

Revisionsdatum: 2024.02.05 Rev. 02

UK Instruction for use for Medicine Cup Model 1010

	Manufacturer	Kellopo A/S Vihelmsborgvej 2 7700 Thisted Denmark	info@kellopo.dk www.kellopo.dk
	Swiss Authorized Representative	The Tao of Excellence GmbH Vorstadt 26 8200 Schaffhausen Switzerland	Kellopo@taoexcellence.ch www.taexcellence.ch
	Importer	Wiegand AG Schlosserstrasse 5 8180 Bülach Switzerland	wiegand@wiegand.ch www.wiegand.ch
	Plastic production	This is the brand of the production. This indicates that it is the plastic production at Kellopo.	
	Production date	The date when the production of the cups has taken place.	
	Shelflife	The medicine cups have a shelf life of 24 months from manufacturing date. The expiry date is printed on the labels, on the medicine cups and directly on the cups. They may be re-used: 1. until the expiry date, or 2. until, the date is no longer visible on the cup whichever is earlier.	
	Storage conditions	Protect from direct sunlight and moisture. Storage temperature of the cups should be as indicated on the symbol – between 5-35°C. Transportation temperatures may vary for the time of the transport, without having any effect on the quality of the cups. The temperature for washing the cups is also higher (up to 65°C), which will also not have an effect on the quality of the cup.	

General
It is well-known worldwide that medication errors are common. This medicine cup is developed to assist trained health professionals in correct dosage and distribution of medicine.
Any serious incident that occurs in connection with the medicine cup should be reported to Kellopo A/S and the competent authority of the Member State in which the user is established.

Pay attention and avoid interruptions when handling medicine (right patient, right drug, right route, right time, right dose, right form of administration within a given route).

Product description
The Medicine Cup is intended for dosage and distribution of non-potent medicine in pill, powder, or liquid form by healthcare professionals.
The cup has volume marks moulded into the plastic, presenting a 25 ml scale and can be washed and re-used by inspecting the presence of date on the cup. If the date is no longer visible, the cup can no longer be used. If the date on the cup has been passed, the cup can also no longer be used.

Warnings
• Make sure that the cup is clean and dry before using it to administer any kind of medicine.
• Ensure the right work surroundings and conditions when dosing liquid medicine (level surface and proper light).
• Check if the printed scale is still visible before the usage.

Measuring scale
The measuring scale is marked from 5ml up to 25ml, in steps of 5ml with an accuracy of +/- 1 ml. To make sure that the scale shows you the proper content, you should only dispense in medicine rooms with good lighting.

Reutilisation
The cups can be washed in a normal washing machine with temperature up to 65°C using normal detergent. The cups can be reutilised until the printed date on the cup is no longer visible, or until the expiry date printed on the cup, whichever happens earlier.
Only cups with a visible date print can be re-used.
Ensure the cups are clean and dry before re-use.

Disposal of the medicine cups
The cups are made of 100% Polypropylene and can be disposed according to the local waste instructions.

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