Packingsbeilage der mylife™ ControlGDH Kontrolllösung

DE

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mylife™ ControlGDH Control Solution Insert

EN

Intended use: With the mylife™ ControlGDH control solution you can perform a quality control check of the mylife™ Unio™, mylife™ Unio™ Cara or the mylife™ Unio™ Neva blood glucose meter and mylife™ Unio™ blood glucose test strips. If the test result is in line with the control solution range the mylife™ Unio™ blood glucose meter and mylife™ Unio™ blood glucose test strips is working correctly. The control solution range is printed on the test strip vial.

The mylife™ ControlGDH control solution is intended for *in vitro* diagnostic (for self-testing outside the body) use only.

Summary and explanation: Control solutions are red liquids containing different levels of glucose (normal, high) that react with the test strips of the BGMS. When used following the same procedure as whole blood, the control solution should give results within the expected range printed on the test strip vial. Control solution should be used whenever you want to check if your BGMS is working correctly. The use of control solution is also recommended for practice testing and to check correct procedure.

Chemical composition: The mylife™ ControlGDH control solutions are aqueous solutions which contain the following ingredients:

	mylife™ ControlGDH-Normal	mylife™ ControlGDH-High
Water	96.822%	97.020%
d-Glucose	0.098%	0.300%
Viscosity enhancing agent	1.000%	1.200%
Inorganic salts/Buffers	1.900%	1.300%
Dye	0.080%	0.080%
Preservative	0.100%	0.100%

* The glucose concentration levels were calibrated using a YSI 2300 analyzer traceable to NIST Standard SRM 917c.

Precautions: The control solution is not intended for human consumption. Do not dilute the control solution with any liquid. Keep the control solution away from children. If the control solution or the vial cap has been swallowed contact your physician immediately.

Storage and handling

- Store your unopened control solution at temperature between 2°C and 30°C (36°F and 86°F).
- Shake the control solution vial before use.
- When you open a new bottle of control solution, please write the opening date on the label. Control solution is good for 3 months after opening the bottle, or until the expiration date printed on the label (whichever comes first).
- Do not freeze.
- Check the expiration date and opening date before use of the control solution.
- Replace the vial cap of the control solution and close tightly right away after use.

Procedure to perform a quality mylife™ ControlGDH control test with the control solution:

1. Select "Control test" in the menu of the meter to perform a test with control solution.
2. Take one test strip from the test strip vial and recap it immediately.
3. Apply the test strip to test strip port of meter.
4. Wait for two seconds until the symbol "Recognising" disappears.
5. Before opening the cap of the control solution, shake the bottle well. Open the bottle and place the cap upright on the table. Drip a drop of control solution on the top of the cap.
6. Gently touch the sample entry of the test strip onto the drop of control solution on the top of the cap.
7. When you hear a beep (if volume is switched on) please wait for the test result of control solution test appearing in five seconds.
8. Clean the top of the cap and replace the cap on the control solution bottle.

Understanding control test results

The results of the quality control test should be within the specified control solution range. This range is printed on the test strip vial label. If the results are within the specified ranges, your meter and test strip work correctly.

Example of control solution range

Control solution range	mg/dL	mmol/L
Normal	83 - 113	4.6 - 6.3
High	255 - 345	14.2 - 19.2

Attention: The values of the control solution mentioned in this leaflet are only for demonstration purpose. Please always refer to the values printed on the label of your test strip vial.

If control solution test results are out of the control range, your BGMS may not be working correctly. Repeat the quality control test considering these instructions and the user manual of the BGMS.

Do not use the BGMS to test your blood glucose if the control solution test results are still outside the range printed on the test strip vial.

The reason your control results could be out of the range:

- Your control solution is expired or has been opened more than 3 months ago.
- Your control solution has been diluted.
- Your test strip is expired. Check the expiration date on the vial of test strips.
- The used test strip is damaged.
- Your control solution or test strip vial was not properly closed or left open for a long time.
- You did not perform the test procedure correctly.
- Malfunction of the meter.
- Your control solution test has been performed out of normal temperature range (< 6°C and > 44°C resp. < 43°F and > 111°F).

Customer service: Please review all the instructions to make sure you are performing the steps correctly. **If you have any questions / problems with mylife™ products, please contact your local distributor's customer service.**

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Emergo Europe B.V.