



Item number 1999915559



## Underpad, ABENA Washable Underpads, 90x85cm, light blue, polyester/rayon/TPU

- ✓ Washable absorbent underpad
- ✓ Superior user comfort
- ✓ OEKO-TEX certified
- ✓ Breathable and liquid safe material
- ✓ Peace of mind and protection around the clock
- ✓ Soft, smooth material



STANDARD 100

### Product Description

ABENA Washable Underpads are primarily used as protective covers for beds and mattresses, providing user peace of mind in any situation. The underpads are made of a washable material that can withstand laundering at temperatures of up to 90 °C. This makes the sheet both a practical and discreet solution for users who want extra protection for incontinence problems. The underpads have a liquid safe barrier layer to prevent leakage. The material is also breathable and the surface is soft and smooth for optimum user comfort. The underpads are designed with rounded corners to reduce wear and tear from regular laundering. ABENA Washable is available in two special versions: with either wings or handles. Our underpads with wings are particularly suitable for active and restless users, as the wings at the sides ensure that the underpad stays put. Alternatively, the other variant with handles at the sides is highly practical in situations where aided repositioning is needed, as the handles make it easy to turn or move the user in their bed.



### Specifications

Base name	Underpad
Brand	ABENA
Sub-brand	Washable Underpads
Sub-color	Light blue
Single or multiple use	Multiple use
Material	Polyester, rayon, TPU
Length/depth	90 cm
Width	85 cm
Absorption Rothwell	2500 ml
Certifications	CE. OEKO-TEX. MD.
CE Class (Medical Devices)	Class I
Product or test standards	ISO 11948-1
Directives, regulations and acts	MDR (EU) 2017/745
Because of the products composition it has an expiry date from the production date on:	5 years
Storage Instructions	Store dry, room temperature and no direct sunlight.
Product Disposal Instructions	Can be disposed of with normal household waste sorted according to local regulations. If contamination is suspected, the product must be discarded as clinical waste.
Packaging Disposal Instructions	Can be recycled or incinerated.

### Instructions for use/application

Ensure that the underpad is placed blue side down on the mattress. Ensure that the underpad is fully smoothed out. For optimum protection, we

recommended washing the product before use. Can be washed up to 100 times.

### Packaging data

Unit	Contains	Length	Width	Height	EAN
cII	30 pcs	48.7 cm	32.4 cm	54.6 cm	5713571007068



Regulation (EU) 2017/745 is a regulation of the European Union on the clinical investigation and sale of medical devices for human use. It repeals Directive 93/42/EEC, which concerns medical devices, and Directive 90/385/EEC, which concerns active implantable medical devices, on 26 May 2021.



The CE mark guarantees that a product is safe to use and complies with all safety precautions. CE stands for Conformité Européenne (European Conformity) and is mainly found on electronic equipment, safety equipment, construction products, and medical devices.



STANDARD 100

OEKO-TEX, also known as "Confidence in Textiles", is a health and environmental label used worldwide. The certificate ensures that textiles do not contain substances that are harmful to health or the environment.



The product is fully breathable for the comfort of the user; also during long-time use. Helps maintain skin integrity.

