



Device type Variations	Stethoscope with single-duct tube, membrane and funnel side, switchable UMDNS code: 13-755 Medical device according to class 1 Stethoscope with dual bell in decorative carton with instructions for use Variations: Aluminium or brass bell Sizes (bell): Ø 47 mm or Ø 33 mm, 9 different colours (tube), glossy or matt
Dimensions	Overall length Adult/child: 730 mm/716 mm Tube length: 550 mm Earpiece length/Ø: 165 mm/5 mm
Bell dimensions	Adult: 47 mm/44 mm/30 mm/22 mm Child: 33 mm/30 mm/30 mm/22 mm (Ø bell/Ø membrane/Ø funnel/height)
Material	Bell: Brass, polished and chrome-plated, or Aluminium, anodised (adult model only) Earpiece: Brass, chrome-plated Yoke spring: Spring steel, chrome-plated Ear buds: TPE/TPS, hard, colour: black Funnel protection ring: NBR Membrane: boPET (Hostaphan RN 250), 0.25 mm thick Moulded tube: Soft PVC, latex-free, resistant against common cleaners and disinfectants
Weight	Adult (brass/aluminium): 164 g/118 g Child (brass): 138 g
Shipping data	Gross weight, packaging & device Adult (brass/aluminium): 264 g/218 g Child (brass): 238 g Box (L x W x H): 280 mm x 98 mm x 37 mm Volume packaging/device: 1015.3 cm ³
Guarantee	2 years
Manufacturer	ERKA. Kallmeyer Medizintechnik GmbH & Co. KG Im Farchet 15 D 83646 Bad Tölz Phone: +49 (0)8041 8009-0 Fax: +49 (0)8041 8009-39 E-mail: info@erka.org
Country of origin	Germany