

Article no. 60091
Category Bib with neck-tie fastening
Product Bib 37 x 70 cm
End use Protection against liquids and food



1. Ideal and hygienic as protection against liquids. Effectively protects the patient.
2. With a welded and reversible pocket, to ensure extra strength and security in case of spills.
3. Easy applied.
4. Reinforced neck-tie, which ensures extra strength and easy application.
5. For a secure fit and to keep well in place, suitable for “fidgety” patients.
6. Made of one layer of tissue and one layer of PE-film.
7. High-quality tissue from 100% virgin fibers, bleached. With high absorption abilities and strength.
8. Strong, resilient and liquid safe white PE-film.
9. The plastic packaging has a perforated opening, which can be re-sealed. Allows a hygienic and practical storage and handling.
10. The packaging is made of white PE.
11. CE-marked product. Fulfils the requirements of (EU) MDR 2017/745.
12. Packed in cartons on pallets. The pallets are stretched wrapped with transparent film.
13. The inner packaging and stretch film could be incinerated during standard conditions, or recycled as standard plastic.
14. The cartons could be incinerated or recycled under standard conditions as cardboard.
15. Produced in Sweden.

Size (cm)	Material	Units/inner package Inner package/transp. unit	Transport units/pallet Height of pallet (m)	Gross/Net weight carton (kg)	Measure of carton LxWxH (cm)	Pallet (EUR) size (mm)
37 x 70	Tissue + PE-film	100 6	54 2,355	5,7 6,3	39 39 24,5	800 x 1200

*Finess Hygiene AB is certified in accordance with ISO 9001:2015, ISO 14001:2015 and ISO 13485:2016.
Our CE-marked products fulfils the requirements of (EU) MDR 2017/745.*

2021-05-26 ver 02