



Technical Data (PD)

Lid for Medicine Cup, Snap Closure				
REF	06-114-0*05			

1.	Description / Specification			
1.01	Description	Lid for Medicine Cup REF 06-103-010x, snap closure		
1.02	Dimension	See Customer Drawing		
1.03	Volume	N/A		
1.04	Material / Resin	Polyethylene (LDPE)		
	Grade of Resin	MAT130		
	Properties of Resin	See www.boettger.com		
1.05	Colour	White		
	Grade of Pigment	PIG125		
1.06	Hygiene Level	Standard hygienic		
1.07	Quality Control	Testing of attributive and variable characteristics in accordance with the valid specification		
1.08	Usage / Purpose	Snap closure of medicine cup		
1.09	Conformities / Regulations	See Classification		
1.10	Product Stewardship	See Product Stewardship of Grades of Resin/Pigment		
1.11	Storage	in closed original packaging, dry, at room temperature		
1.12	Transport	Dry, at room temperature, handle with care, do not throw or push cartons, do not stack pallets		
1.13	Other information	For single use only, recyclable		


















2.	Features			
2.01	Basic features	-		
2.02	Temperature range	-50°C to 60°C		
2.03	Autoclavability	No		
2.04	Centrifugation	N/A		
2.05	Chemical Resistance	See www.boettger.com		
2.06	Expiry	5 years from LOT (YYYY-MM)		
2.07	Other information	-		










3.	Packaging	06-114-0005	06-114-0105	
3.01	Pcs/Bag	750	75	
3.02	Pcs/Ctn	3750	3750	
3.03	Pcs/Pallet	180000	90000	
3.06	LOT / VP	WWYY (week, year of production) / XX (packaging)		
3.07	Other information	See Measures		

4.	Other Information			
4.01	Codes and Numbers	See Codes and Numbers		
4.02	Symbols	See Symbols		
4.03	CH-REP	Swiss AR Services AG, Industriestr. 47, CH-6300 Zug CHRN-AR-20000807		

Revision	Drawn	Approved	Confidential: Information contained in this document or drawing is confidential and proprietary to Paul Boettger GmbH & Co. KG. This document may not be reproduced for any reason without written permission from Paul Boettger GmbH & Co. KG. All rights of design, invention and copyright are reserved.
16	04.03.2024 QMB	04.03.2024 QMB	

Technical Data (PD)

5.	Label of Sales Unit (Example) according to MDR 2017/745/EC															
	<div style="display: flex; justify-content: space-between; align-items: flex-start;"> <div style="width: 60%;"> <p>Paul Boettger GmbH & Co. KG Miesleuthenweg 35 - 45 94249 Bodenmais/Germany www.boettger.com</p> </div> <div style="width: 35%; text-align: right;">  boettger </div> </div> <hr/> <div style="display: flex; justify-content: space-between; align-items: center;"> 3750 St/pcs REF 06-114-0005 </div> <hr/> <p>Deckel Medizinbecher PE (weiß) Lid Medicine Cup PE (white)</p> <div style="text-align: right; font-size: small;">V240201</div> <hr/> <div style="display: flex; justify-content: space-between; align-items: center;"> <div style="width: 30%;"> <table border="1" style="font-size: x-small;"> <tr> <td style="text-align: center;">MD</td> <td rowspan="2" style="text-align: center;">  </td> <td rowspan="2" style="text-align: center;">  </td> <td rowspan="2" style="text-align: center;"> Paul Boettger GmbH & Co KG 94249 Bodenmais Germany </td> <td rowspan="2" style="text-align: center;">  </td> <td rowspan="2" style="text-align: center;"> 2099-99 </td> <td rowspan="2" style="text-align: center;">  </td> </tr> <tr> <td style="text-align: center;">CH REP</td> <td style="text-align: center;"> Swiss AR Services AG, Industriestrasse 47, CH-6300 Zug </td> <td style="text-align: center;"> Made in Germany </td> </tr> </table> </div> <div style="width: 35%; text-align: center;">  </div> </div> <hr/> <div style="display: flex; justify-content: space-between; align-items: center;"> <div style="width: 30%;"> <table border="1" style="font-size: x-small;"> <tr> <td style="text-align: center;">UDI</td> <td style="text-align: center;">  (01) 0 4042836 00714 9 (17) 999900 (10) 0000 </td> </tr> </table> </div> <div style="width: 35%; text-align: right;"> <table border="1" style="font-size: x-small;"> <tr> <td style="text-align: center;">LOT</td> </tr> <tr> <td style="text-align: center;">VP</td> </tr> <tr> <td style="text-align: center;">Ctn 1</td> </tr> </table> </div> </div>	MD			Paul Boettger GmbH & Co KG 94249 Bodenmais Germany		2099-99		CH REP	Swiss AR Services AG, Industriestrasse 47, CH-6300 Zug	Made in Germany	UDI	 (01) 0 4042836 00714 9 (17) 999900 (10) 0000	LOT	VP	Ctn 1
MD									Paul Boettger GmbH & Co KG 94249 Bodenmais Germany		2099-99					
CH REP			Swiss AR Services AG, Industriestrasse 47, CH-6300 Zug	Made in Germany												
UDI	 (01) 0 4042836 00714 9 (17) 999900 (10) 0000															
LOT																
VP																
Ctn 1																

6.	Printing on Subpackaging (Example)									
	<div style="display: flex; justify-content: space-between; align-items: center;"> <div style="width: 30%;"> <table border="1" style="font-size: x-small;"> <tr> <td style="text-align: center;">REF</td> <td style="text-align: center;">06-114-0*05</td> <td style="text-align: center;">  </td> </tr> <tr> <td style="text-align: center;">  </td> <td style="text-align: center;">2021-10-19</td> <td></td> </tr> <tr> <td style="text-align: center;">  </td> <td style="text-align: center;">Paul Boettger GmbH & Co. KG</td> <td></td> </tr> </table> </div> </div>	REF	06-114-0*05			2021-10-19			Paul Boettger GmbH & Co. KG	
REF	06-114-0*05									
	2021-10-19									
	Paul Boettger GmbH & Co. KG									

* Placeholder for packaging alternatives

All informations without guarantee. Above mentioned products are manufactured according to company specifications and our certified Quality Management System corresponding to DIN EN ISO 13485. The quality check is carried out on a sample basis; an individual examination is not possible. Deficiencies parts may occur. We accept no liability for consequential damages.



Revision	Drawn	Approved	Confidential: Information contained in this document or drawing is confidential and proprietary to Paul Boettger GmbH & Co. KG. This document may not be reproduced for any reason without written permission from Paul Boettger GmbH & Co. KG. All rights of design, invention and copyright are reserved.
16	04.03.2024 QMB	04.03.2024 QMB	