



boettger

Technical Data (TD)

Medicine Cup				
REF	06-103-010*			

1.	Description / Specification			
1.01	Description	Universal Container 30		
1.02	Dimension	See Customer Drawing		
1.03	Volume	No Measurement Functionality		
1.04	Material / Resin	Polypropylene (PP)		
	Grade of Resin	MAT111		
	Properties of Resin	See www.boettger.com		
1.05	REF	06-103-0100	06-103-0101	06-065-0102
	Colour	Natural	Blue	Red
	Grade of Pigment		PIG101	PIG102
			06-103-0103	06-065-0104
			Yellow	Green
	PIG103	PIG104		
1.06	Hygiene Level	Standard hygienic		
1.07	Quality Control	Testing of attributive and variable characteristics in accordance with the valid specification		
1.08	Usage / Purpose	for collecting / distributing tablets		
1.09	Conformities / Regulations	See Classification		
1.10	Product Stewardship	See Product Stewardship of Grades of Resin/Pigment		
1.11	Storage	in closed original packaging, dry, at room temperature		
1.12	Transport	Dry, at room temperature, handle with care, do not throw or push cartons, do not stack pallets		
1.13	Other information	For single use only, recyclable		

2.	Features			
2.01	Basic features	-		
2.02	Temperature range	-50°C to 60°C		
2.03	Autoclavability	No		
2.04	Centrifugation	N/A		
2.05	Chemical Resistance	See www.boettger.com		
2.06	Expiry	5 years from LOT (YYYY-MM)		
2.07	Other information	-		

3.	Packaging	06-103-010*		
3.01	Pcs/Bag	75		
3.02	Pcs/Ctn	3750		
3.03	Pcs/Pallet	180000		
3.04	LOT / VP	WWYY (week, year of production) / XX (packaging)		
3.05	Other information	See Measures		


Revision	Drawn	Approved	Confidential:
14	04.03.2024	04.03.2024	Information contained in this document or drawing is confidential and proprietary to Paul Boettger GmbH & Co. KG. This document may not be reproduced for any reason without written permission from Paul Boettger GmbH & Co. KG. All rights of design, invention and copyright are reserved.
	QMB	QMB	




boettger

Technical Data (TD)

4.	Other Information	
4.01	Codes and Numbers	See Codes and Numbers
4.02	Symbols	See Symbols
4.03	CH-REP	Swiss AR Services AG, Industriestr. 47, CH-6300 Zug CHRN-AR-20000807

5.	Label of Sales Unit (Example) according to MDR 2017/745/EC	
	<p>Paul Boettger GmbH & Co. KG Miesleuthenweg 35 - 45 94249 Bodenmais/Germany www.boettger.com</p>  <hr/> <p>3750 St/pcs REF 06-103-0100</p> <hr/> <p>Medizinbecher PP (natur), Medicine Cup PP (natural) zum Sammeln/Verteilen von Tabletten</p> <p style="text-align: right;">V240201</p> <hr/> <p>    Paul Boettger GmbH & Co KG 94249 Bodenmais Germany  2099-99   </p> <hr/> <p>  Swiss AR Services AG, Industriestrasse 47, CH-6300 Zug  Made in Germany </p> <hr/> <p>    </p> <p style="text-align: center;">(01) 0 4042836 00199 4 (17) 999900 (10) 0000</p> <p style="text-align: right;"> VP Ctn 1 </p>	

6.	Printing on Subpackaging (Example)	
	<p>  06-103-010*  </p> <p>  2021-10-19 </p> <p>  Paul Boettger GmbH & Co. KG </p>	

* Placeholder for colour alternatives

All informations without guarantee. Above mentioned products are manufactured according to company specifications and our certified Quality Management System corresponding to DIN EN ISO 13485. The quality check is carried out on a sample basis; an individual examination is not possible. Deficiencies parts may occur. We accept no liability for consequential damages.



Revision	Drawn	Approved	Confidential: Information contained in this document or drawing is confidential and proprietary to Paul Boettger GmbH & Co. KG. This document may not be reproduced for any reason without written permission from Paul Boettger GmbH & Co. KG. All rights of design, invention and copyright are reserved.
14	04.03.2024 QMB	04.03.2024 QMB	