



2365 Camino Vida Roble
Carlsbad, CA 92011
United States

Date: 17-March-2017

Reference: 1616, Adult Nasal Oxygen Cannula with 3-channel tubing.

Data Sheet

| Part Description | Material |
|----------------------|--------------------------|
| Facepiece (Adult) | Polyvinylchloride |
| Tubing, 19" | Polyvinylchloride |
| Bolo | Low Density Polyethylene |
| Connector, "Y" | Polyvinylchloride |
| Tubing, 3-channel | Polyvinylchloride |
| Connector, End Piece | Polyvinylchloride |

Salter Labs manufactures product from medical grade materials in compliance with Good Manufacturing Practices (GMP's) as listed in 21 C.F.R. (U.S. code of Federal Regulations)

I certify that the above information is accurate to the best of my abilities.

Please do not hesitate to contact me if I can be of any further assistance. My email address is mcollins@salterlabs.com and my phone number is 760-795-7102.

Peggy Collins
Regulatory Affairs
Salter Labs