



2365 Camino Vida Roble
Carlsbad, CA 92011 United
States

Date: 5-September-2018

Reference: 16SOFT, Nasal Oxygen Cannula with 3-channel tubing.

Data Sheet

Part Description	Material
Facepiece (Adult)	Polyvinylchloride (PVC)
Tubing, 19"	Polyvinylchloride (PVC)
Bolo	Low Density Polyethylene
Connector, "Y"	Polyvinylchloride (PVC)
Tubing, 3-channel	Polyvinylchloride (PVC)
Connector, End Piece	Polyvinylchloride (PVC)

Nasal Cannulas with PVC oxygen tubing have been in use for more than 40 years and are an industry standard. The materials used in the Salter Labs cannulas have been tested and found to meet the requirements of ISO 10993 for biocompatibility (irritation, sensitization, and cytotoxicity).

Salter Labs manufactures product from medical grade materials in compliance with Good Manufacturing Practices (GMP's) as listed in 21 C.F.R. (U.S. code of Federal Regulations).

Please do not hesitate to contact me if I can be of any further assistance. My email address is mcollins@Salterlabs.com and my phone number is 760-795-7102.

Peggy Collins
Regulatory Affairs