



2365 Camino Vida Roble
Carlsbad, CA 92011 United
States

Date: 17-March-2017

Reference: 1053, Conventional Nasal Oxygen Cannula with 3 channel tubing.

Data Sheet

Part Description	Material
Facepiece (Adult)	Polyvinylchloride
Tubing, . 0.4826 m	Polyvinylchloride
Bolo	Low Density Polyethylene
Connector, "Y"	Polyvinylchloride
Tubing, 3-channel	Polyvinylchloride
Connector, End Piece	Polyvinylchloride

Salter Labs manufactures product from medical grade materials in compliance with Good Manufacturing Practices (GMP's) as listed in 21 C.F.R. (U.S. code of Federal Regulations)

I certify that the above information is accurate to the best of my abilities.

Peggy Collins
Regulatory Affairs

Peggy Collins | Senior Regulatory Affairs Specialist | **Salter Labs**
2365 Camino Vida Roble
Carlsbad, CA 92011
[760-795-7102](tel:760-795-7102)
www.salterlabs.com
Design Innovations for Respiratory Care

www.salterlabs.com

Telephone: 760-795-7100 Fax: 760-683-6797