



2365 Camino Vida Robles  
Carlsbad, CA 9201  
United States

**Date:** 16-November-2016

**Reference No.:** 8130, Adult Oxygen Mask, High Concentration, Non-Rebreathing with 7 ft. 3-channel tubing, and safety vent

Part Description	Material
Mask, Assembly	
Mask, Medium Concentration (Adult)	Polyvinylchloride
Safety Vent	Polyvinylchloride
Nose Bar	Aluminum
Elastic Strap	Elastic & Cloth
Reservoir Bag Assembly	
Tubing, 3-channel, 7 ft	Polyvinylchloride
Connector, End Piece	Polyvinylchloride
Rebreather	Polyvinylchloride
Tape, Bag	Polyvinylchloride
Grommet, Oxygen	Polyvinylchloride (Rigid)

Salter Labs manufactures product from medical grade materials in compliance with Good Manufacturing Practices (GMP's) as listed in 21 C.F.R. (U.S. code of Federal Regulations)

I certify that the above information is accurate to the best of my abilities.

Please do not hesitate to contact me if I can be of any further assistance. My email address is [mcollins@salterlabs.com](mailto:mcollins@salterlabs.com) and my phone number is 760-795-7102.

Regulatory Affairs