



2365 Camino Vida Roble
Carlsbad, CA 92011
United States

Date: 17-March-2017

Reference: SO-407, Adult Oxygen Mask Medium Concentration, Supply Tube 7' (2.1 m).

Data Sheet

Part Description	Material
Mask, Assembly	
Mask, MediumConcentration (Adult)	Polyvinylchloride
Safety Vent	Polyvinylchloride
Nose Bar	Aluminum
Elastic Strap	Elastic & Cloth

Salter Labs manufactures product from medical grade materials in compliance with Good Manufacturing Practices (GMP's) as listed in 21 C.F.R. (U.S. code of Federal Regulations)

I certify that the above information is accurate to the best of my abilities.

Please do not hesitate to contact me if I can be of any further assistance. My email address is mcollins@salterlabs.com and my phone number is 760-795-7102.

Regulatory Affairs