

# ALCO-PREP® Large

## Cleansing of the intact, healthy skin



### HOS-610HD



#### Description

- Prep pad fully saturated with 70% isopropyl alcohol U.S.P.
- Hermetically sealed in a laminated aluminium foil pouch.
- For single and external use only.

#### Application

- Cleansing of the intact, healthy skin

#### Commercial Data

- Item ID HOS-610HD
- Type ALCO-PREP® Large - Cleansing pads for external use
- Basic UDI-DI 542501868ALCOPREPLARGE0AX
- CNK code (Belgium) 3466-224

#### Technical Data

- Dimensions pad 34 x 68 mm (unfolded)
- Material Soft folded, pin doted, non-woven prep pad fully saturated with 70% isopropyl alcohol U.S.P.
- Validity 5 years
- Marking CE-marked, Class 1, non-sterile  
Lot number and expiration date printed on each pouch, box and case
- Material Safety Data Sheet Available upon request
- Certificate of Analysis Available with every batch produced upon request
- Quality Management Hospidex NV / Hospidex France SAS are ISO EN 13485:2016 certified

#### Packaging

	<b>Box – 100 pouches</b>	<b>Case - 20 boxes</b>	<b>Pallet - 100 cases</b>
• Dimensions	10,4 x 6,8 x 6,0 cm	54 x 15 x 14 cm	120 x 80 x 155 cm
• Weight	0,120 kg	2,8 kg	295 kg