

ERKAphone Duo

**Straightforward. Universal.**

The ERKAphon product family offers an extensive selection of models for the various areas of application of a stethoscope.

Available for adults and children, the ERKAphon Duo features a switchable double bell. The open funnel side is designed especially for low frequencies, the membrane side for the entire frequency range. The interplay of high-end components in the ERKAphon series ensures the best possible acoustics.



Since the ERKA standard membranes are installed with convex basic stress, they guarantee the best hearing. A ring as protection against cold is available for the membrane side as an option.



The bell – available in brass/chrome-plated brass with high-gloss finish/chrome-plated with matt finish or anodised aluminium – is quick and easy to clean.



The tube is available in nine colour variations.



To ensure that signals captured by the bell are transmitted to the ears without interference, the spring was not integrated in the tube duct.



The ERKAphon Duo is also available in a special version for teaching applications.