



**Salvaclip Standard perineal clip belt** to help positioning of seated patients. Especially useful for people with stability problems due to a functional disability. Perineal support prevents the patient from moving downward.

Side release plastic buckles allows a quick and easy application and can easily adjust the product depending on the size of the patient. It is the best choice for use in patients at risk of slipping and falling off the chair, those who adopt incorrect postures or have uncontrolled movements.


Made of polyester with water repellent coating and comfortable polyurethane foam. Reinforced polyester straps.

## PERINEAL CLIP BELT


### FEATURES



Side release plastic buckles



Anti stain liquid repellent



Reinforced polyester straps



Suitable for wheelchair and armchair

### SPECIFICATIONS

PRODUCT FAMILY: SALVACLIP Standard

PRODUCT NAME: Perineal clip belt

#### SIZE CHART

Ref.	Size	Measures (cm)
		Contour Min./Max.
SCP4330	M	90 / 165
SCP4330L	L	115 / 200

#### COLOUR

- Ext.Blue - Int.Black
- Dark blue straps

#### COMPOSITION


**Belt:** Ext: Polyester - Int: Polyurethane  
**Straps:** Polyester  
**Buckles:** Polypropilene (PP)

#### LAUNDRY INSTRUCTIONS



- Machine wash 40°C
- Do not bleach
- Tumble dry low
- Do not iron
- Do not dry clean

#### PACKAGING



- 1 - Salvaclip Standard - Perineal clip belt
- 1 - Safety and instructions manual