

Product Identification

Product family	Sterile Dressing sets
Reference	1040042
Name	Woundcare pack No. 017 *S
Sterilization	Ethylene Oxide
Brand	Sterisets®



Material

Code	Designation	Qty	Composition
5015010	Ball NW 13 x 15 cm, Ø 35 mm, 30 g/m ² , with ring	5	Non woven 30 g/m ²
5021605	Drape 2 ply PE 25 gr/m ² NW 30 gr/m ² 40 x 45 cm	1	PE/NW 55g/m ² (polyethylene 27 g/m ² + non woven 28 g/m ²)
8100110	Forceps Burger 12,6 cm plastic green	1	Polypropylene 30% glass filled
8100114	Forceps Burger 12,6 cm plastic blue	1	Polypropylene 30% glass filled

Packaging

Type	Sterile single pack
Package weight	34 g
Package dimensions	165 x 125 x 20 mm
Package UDI-DI (GTIN)	05608120026129
Quantity per inner box	Not applicable
Inner box UDI-DI (GTIN)	Not applicable
Quantity transport carton	40 pcs
Transport carton UDI-DI (GTIN)	15608120026126



Hazards Identification

Latex	No
Medicinal substances	No
Animal tissue	No
Toxic or dangerous materials	No

Shelf life

5 years from the date of manufacture, unless any component has a shorter shelf-life.

Storage conditions

Keep away from sunlight and keep dry.

Environmental and disposal considerations

Not regarded as dangerous for the environmental. Dispose of material according to local, regional and national regulations.

Compliance Regulations & Standards

Regulation (EU) 2017/745 (as applicable), Directive 93/42/EEC, ISO 13485, ISO 14971, ISO 15223-1, ISO 10993-1, ISO 11737, ISO 11135.

