

Product Identification	
Product family	Sterile Dressing sets
Reference	1040044
Name	Woundcare pack No. 018 *S
Sterilization	Ethylene Oxide
Brand	Sterisets®



Material			
Code	Designation	Qty	Composition
5012035	Swab 7,5 x 7,5 cm, 30 gr/4 ply NW - per 5	5	Non woven 30 g/m ²
5015010	Ball NW 13 x 15 cm, Ø 35 mm, 30 g/m ² , with ring	5	Non woven 30 g/m ²
5021605	Drape 2 ply PE 25 gr/m ² NW 30 gr/m ² 40 x 45 cm	1	PE/NW
8100110	Forceps Burger 12,6 cm plastic green	1	Polypropylene, 30% glass filled
8100114	Forceps Burger 12,6 cm plastic blue	1	Polypropylene, 30% glass filled

Packaging	
Type	Sterile single pack
Package weight	38 g
Package dimensions	165 x 125 x 25 mm
Package UDI-DI (GTIN)	05608120025474
Quantity per inner box	Not applicable
Inner box UDI-DI (GTIN)	Not applicable
Quantity transport carton	40 pcs
Transport carton UDI-DI (GTIN)	15608120025471

Hazards Identification

Latex	No
Medicinal substances	No
Animal tissue	No
Toxic or dangerous materials	No

Shelf life

5 years from the date of manufacture, unless any component has a shorter shelf-life.

Storage conditions

Keep away from sunlight and keep dry.

Environmental and disposal considerations

Not regarded as dangerous for the environmental. Dispose of material according to local, regional and national regulations.

Compliance Regulations & Standards

Regulation (EU) 2017/745 (as applicable), Directive 93/42/EEC, ISO 13485, ISO 14971, ISO 15223-1, ISO 10993-1, ISO 11737, ISO 11135.

