


ULCER X Masstabellen



Umfang in cm	x-small	small	small PLUS	medium	medium PLUS	large	large PLUS	x-large	x-large PLUS
cD	26-32	28-34	32-38	31-37	36-42	34-40	40-46	38-43	44-49
cC	31-36	31-36	35-40	34-39	39-44	37-42	43-48	40-45	47-52
cB'	24-29	24-29	24-29	29-32	29-32	32-36	32-36	35-39	35-39
cB	17-19	19-23	19-23	23-26	23-26	26-29	26-29	29-32	29-32
cY	26-32	26-32	28-32	32-35	32-35	35-38	35-38	38-41	38-41
cA	17-19	19-23	19-23	23-26	23-26	26-29	26-29	29-32	29-32

c = Umfang
ℓ = Länge

Länge in cm	
kurz	33-37
lang	38-44